# ARTIKULATOR CA 3.0

INSTRUCTIONS FOR USE

---

GEBRAUCHSINFORMATION

ANWENDUNG  
 APPLICATION  
 UTILISATION  
 USO  
 APLICACIÓN  
 GEBRUIK  
 ПРИМЕНЕНИЕ

1.0

Direktes Eingipsen dank Split-Cast am Artikulator  
 Direct fixing with plaster thanks to the Split-Cast system  
 Plâtre direct avec socle Split-Cast sur l'articulateur  
 Cessatura diretta grazie allo Split-Cast collegato all'articulatore  
 Enyesado directo gracias al sistema Split-Cast en el articulador  
 Direct ingipsen dankzij Split-Cast op de articulator  
 Прямая гипсовка благодаря разборному цоколию (split-cast) на артикуляторе

2.0

Modellplatten Kunststoff, Modellplatten Metall  
 Plastic mounting plates, Metal mounting plates  
 Plaques modèles : en matière synthétique, en métal  
 Piastre per modelli: in plastica, in metallo  
 Placas de modelo de plástico, Placas de modelo metálicas  
 Modelplaten in kunststof, Metalen modelplaten  
 Пластмассовые монтажные пластины, металлические монтажные пластины

6.0

Protrusionsbewegung  
 Protrusive movement  
 Mouvement de protrusion  
 Movimento di protrusione  
 Movimiento protrusivo  
 Protrusiebeweging  
 Протрузионное движение

7.0

Lateral- und Bennettbewegung  
 Lateral and Bennett movement  
 Mouvement latéral et de Bennett  
 Movimento di lateralità e di Bennett  
 Movimiento de lateralidad y de Bennett  
 Laterale beweging en beweging van Bennett  
 Латеральное движение и движение Беннета

3.0

Geräteinstellung Mittelwert 30°  
Initial setting up, average value of 30°  
Valeur moyenne de réglage de l'appareil 30°  
Regolazione strumento valore medio 30°  
Ajuste del aparato con valor medio 30°  
Instelling apparaat gemiddelde waarde 30°  
Настройка прибора, среднее значение 30°

4.0

Zentrixfixierung geschlossen/open  
Centric lock closed/open  
Fixation centrée fermée/ouverte  
Fissaggio di centrica chiuso/aperto  
Fijación céntrica cerrada/abierta  
Centrische snelsluiting gesloten/open  
Центральный фиксатор закрыт/открыт

5.0

Zentrixfixierung offen. Retrusionsriegel offen  
Centric lock open. Retraction lock open  
Fixation centrée ouverte. Verrouillage de rétrusion ouvert  
Fissaggio di centrica aperto. Blocco di retrusione aperto  
Fijación céntrica abierta. Cierre de retrusión abierto  
Centrische snelsluiting open. Retrusiestaaf open  
Центральный фиксатор открыт. Ретрузионный винт открыт

8.0

Immediate Side Shift (ISS) Bewegung  
Immediate side shift (ISS)  
Déplacement immediate side shift (ISS)  
Movimento immediate side shift (ISS)  
Movimiento immediate side shift (ISS)  
Immediate side shift (ISS)-beweging  
Движение непосредственного бокового сдвига (ISS)

9.0

Retrusionsbewegung  
Retrusive movement  
Voie de rétrusion  
Movimento di retrusione  
Movimiento retrusivo  
Retrusiebeweging  
Ретрузионное движение

10.0

O-Ring Montage  
Fitting the O-ring seal  
Montage du joint torique  
Montaggio anello ad O  
Montaje del anillo toroidal  
Montage van de rubberen ring  
Установка O-образного кольца

ERSATZTEILE  
SPARE PARTS  
PIÈCES DE RECHANGE  
ELENCO PARTICOLARI  
REPUESTOS  
WISSELSTUKKEN  
ЗАПАСНЫЕ ЧАСТИ

# ARTIKULATOR CA 3.0

**DE** GEBRAUCHSINFORMATION

**EN** INSTRUCTIONS FOR USE

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# INSTRUCTIONS FOR USE – ARTICULATOR CA 3.0

---

## Specifications

- > Bonwill's triangle 110 mm
  - > Balkwill angle 25°
  - > Height 150 mm, inside height 100 mm with mounting plates, 115 mm without, depth 160 mm, width 145 mm
  - > Retrusive path 1.5 mm
  - > Immediate side shift (ISS): 0–2.5 mm fixed
  - > Condyle path inclination adjustable 0°–60°
  - > Incisal guide 15°
  - > Weight: 690 g
  - > Material: anodised aluminium
  - > Colour: silver
- 

## 1. Introduction

### 1.1 Preface

Thank you for purchasing a branded product from CANDULOR. This appliance is a technically sophisticated product. Top quality and high precision are the hallmarks of this equipment. If used improperly, parts of the CA 3.0 may be damaged, so please read the instructions for use before using it for the first time. We hope you will enjoy working with the CA 3.0.

### 1.2 Warranty

We provide a 2-year warranty from the date of purchase against material failure of original metal parts. The validity of the warranty depends on proper handling in accordance with the Instructions for Use.

## 2. Safety

This appliance has been developed solely for use in dentistry and must be handled strictly in accordance with the Instructions for Use. The manufacturer cannot be held liable for any damages arising as a result of incorrect handling or of using the appliance for other purposes. The user shall be solely responsible for testing the appliance prior to use in order to verify its suitability for the intended purpose, particularly if that purpose is not explicitly stated in the instructions.

## 3. Product Description

### 3.1 Function description, layout of the appliance

The CA 3.0 is a semi-adjustable articulator designed in accordance with Camper's plane. The mechanical joint takes the form of a guide edge which supports a bi-conical axle. The horizontal

condyle path inclination can be adjusted from 0–60° by rotating a wheel. The appliance incorporates a retrusive movement of 1.5 mm as well as the possibility of an immediate side shift. A novel centric quick-lock facility can be activated for exact, reproducible initial set position of the articulator. The upper and lower sections of the articulator remain connected by the centric rubber bands even when the centric lock is open. Thanks to its well-conceived range of accessories the CA 3.0 can be readily adapted to suit the personal preferences and requirements of the user.

### 3.2 Indication

Suitable for the three-dimensional fixing of models of dental reconstructions. The basic configuration is suitable for simulating the average-value movements of the human temporomandibular joint.

### 3.3 Contraindications

Based on current knowledge, there are no contraindications providing that the appliance is handled properly in accordance with the Instructions for Use.

## 4. Installation and First Use

### 4.1 Using the appliance for the first time

Remove the parts from the box and check for completeness. If any parts are missing or damaged, please contact Customer Service. We recommend that you keep the packaging for transportation purposes. Before using the articulator for the first time, spray with a standard silicone spray. Do not grease!

### 4.2 Setting up the appliance (average value) (3.0)

- > Loosen knurled thumb screw A and set condyle path inclination to an average value of 30°.
- > Close the centric lock. (4.0)
- > Push the vertical guide rod on the vertical guide rod holder into the zero position (millimetre scale) and fasten.
- > Centre the incisal guide on the vertical guide rod
- > Push the long part of the incisal pointer through the top hole and the short part through the bottom hole on the vertical guide rod as far as they will go.

### 4.3 Support rod

The supporting rod is a feature of the ergonomic design of the articulator. Screw the supporting rod

into the thread in the upper section of the CA 3.0. When the supporting rod is fitted the upper section of the articulator will be supported (opening stop).

### 4.4 Mounting in the articulator (1.0 + 2.0)

The following possibilities are available for mounting in the articulator:

- > Direct fixing with plaster thanks to the split-cast system (1.0)
- > Plastic mounting plates, Metal mounting plates (2.0)

## 5. Handling, Operation

### 5.1 Centric position

The CA 3.0 features a precision centric lock which can be activated in 3 possible positions:

#### Centric lock (4.0) closed

In this position the CA 3.0 can only be opened and closed in centric position for occlusion checks to be performed.

#### Centric lock (4.0) open

#### Retrusion lock closed

Protrusive, Bennett, balance movements as well

as immediate side shift (ISS) and all intermediate movements can be physiologically reproduced in this position.

#### Centric lock and retrusion lock open (5.0)

With the lock in this position, the retrusive movement can also be performed.

### 5.2 Protrusive movement (6.0)

If the centric lock is open, the protrusive movement can be performed as shown in the figure. The horizontal condyle path inclination can be adjusted from 0–60° by rotating the wheel.

### 5.3 Lateral and Bennett movement (7.0)

If the centric lock is open, the lateral movements can be performed as shown in the figure. This is done by applying backward pressure with the thumb on one side on the appropriate joint.

### 5.4 Immediate side shift (ISS) (8.0)

If the centric lock is open, the immediate side shift can be performed as shown in the figure. This is done by pressing in the transverse direction on one side with the index finger.

### 5.5 Retrusive movement (9.0)

A retrusive movement can only be performed if the centric lock and the retrusion lock are open. This opens the way for retrusive movement.

## 6. Practical Application of Model Orientation

### 6.1 Dentulous or edentulous case with rubber band (average value)

- > Set up the appliance as described under 4.2.
- > The notches in the side columns and the circular groove on the vertical guide rod can be used to orient the occlusal plane.
- > Place the plaster on the mounting plate and base block, and slowly close the articulator. Next, fit a thin rubber band. Then attach it around the columns, in the notches and beneath the longer part of the incisal pointer.
- > Align the occlusal plane with the rubber band.
- > Align the incisal pointer between the mandibular central incisors or on the wax rim with the respective mark. Ideally, fix the lower model with modelling material and align. Finally, cast the upper model in plaster.

## 7. Maintenance

### 7.1 Cleaning

The CA 3.0 is a precision instrument and should only be steam cleaned. If no steam-jet equipment is available, it should be rinsed off with clean, hot water. All contact with strong acids and solvents and cleaning with ultrasonic equipment are to be avoided. Do not scrape the surface of the appliance with sharp objects.

### 7.2 Maintenance checks (10.0)

How often these checks need to be carried out will depend on the intensity of use and the user's working methods. For this reason, the stated frequency should only be taken as a guide.

*Check* – for cracks and signs of damage

*Part* – O ring seal

*Frequency* – Once a month

### *How to replace*

- > Open centric lock, and open and tighten retrusion lock.
- > With the articulator open and lying on its side, slip the O ring seal over the joint axle from the outside and push over the guide edges using a blunt instrument.
- > Fit O-ring seal. (10.0)

### 7.3 Repairs

The warranty will become null and void if any attempts at repair are performed within the warranty period by anyone other than a qualified CANDULOR service centre. Only original CANDULOR spare parts and accessories may be used.

## 8. Delivery Scope (Art. No. 668314)

- 1 × CA 3.0
- 1 × Incisal guide 15°
- 1 × Incisal pointer
- 2 × Retention discs
- 2 × Plastic mounting plates
- 1 × Supporting rod for opening stop
- 2 × Plastering covers
- 1 × Instructions for use

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**CANDULOR AG**