# Biodent<sup>®</sup> K+B Plus

Instructions for Use  
Gebrauchsinformation  
Mode d'emploi  
Istruzioni per l'uso  
Bruksanvisning

**DENTSPLY**

Fig. 1


Fig. 2

Fig. 3


Fig. 4


Fig. 5


Fig. 6


Fig. 7

Fig. 8

Fig. 9


Fig. 10


Fig. 11

## Applications

- Crown and bridge veneers
- Jacket crowns in Biodent shades
- Attachment of plastic teeth in the flask pack method
- Temporary crowns and bridges
- Corrections
- Repairs
- Telescope and cone works

K+B Plus is not suitable for occlusal surfaces.

## Technical Data

Cured in pressure vessel (Type 1 according to DIN EN ISO 10477:1995).

Self-curing (Type 2 according to DIN EN ISO 10477:1995).

Main components:

Powder: Polymethyl methacrylates  
Liquid: Methyl methacrylates

In accordance with DIN EN ISO 10477:1995

1. Bending strength: 65 MPa
2. E module: 2000 MPa
3. Water uptake: complies with standard specification
4. Water solubility: complies with standard specification
5. Colour stability: complies with standard specification

## Characteristics

### Material

Biodent K+B Plus is made of the dental plastic INPEN, a blend of selected polymers without filler and with a high proportion of cross-linked particles. The proportions of the various ingredients were formulated by computer to produce an ideal composition. The result was "Interpenetrated Polymer Network" = INPEN – a heterogeneous micro-structure.

Advantages: crowns and veneers resistant to abrasion, tough, yet not brittle; no problems with grinding; easy to bring quickly to a very high polish; permanently thick and smooth surfaces; uniform curing temperature; reliable bonding of additional layers.

### Colour

Biodent shades are easily and reliably obtained in the plastic veneer. Each shade has its own dentine opaque and the shading of the bevelled layers is supported by nine different neck opaques. The sheen and natural qualities of the finished Biodent K+B workpieces have a further optical quality of natural teeth added to them: white fluorescence under ultra-violet light!

### Technique

Three different liquids enable the use of any of the following techniques:

#### **Liquid S for the direct build-up method**

Short doughing time (2–4 minutes at 22–30 °C). Excellent plasticity, i.e. very easy to mould and cut (not rubbery). The dough is not tacky and does not cobweb. The materials do not run into each other during the layering but set quickly. No drying out, no formation of skin. Ample working time (45 minutes).

Cure: water bath at 95 °C, at least 15 minutes, pressure 6 bars.

#### **Liquid N for flask method**

Good fluidity, long working time.

Cure: heat up for about 30 minutes, then boil for 30 minutes.

#### **Liquid K for cold cure**

Sets quickly, short hardening time, only in conjunction with Promotor K.

Cure: 8–12 minutes (according to room temperature) or 10 minutes in pressure flask at 40–50 °C.

## OVS Opaque II

Liquid, ready to use, covers perfectly in a very thin film.

Drying time: 15 minutes at room temperature. This can be reduced to 5 minutes by warming.

## Presentation

### Single Items

Dentine-, Enamel powders	20 g, 100 g
Neck-, Transparent-, Tissue tone powders	20 g

### Liquid S

Liquid S	30 ml, 100 ml
Liquid N	30 ml, 100 ml
Liquid K	25 ml, 80 ml

### Accessories

Protective varnish	30 ml
Protective varnish thinner	10 ml
Promotor K	5 g, 15 g, 50 g
Mixing vessel for Promotor K	
Biodent Isolation Fluid H	15 ml
Biodent Isolation Fluid K	30 ml
Mixing palette with glass lid	
Set of brush tips (10 pieces)	
Measuring spoon	

### OVS II Opaque Kit

**For Biodent K+B Plus use exclusively the Biodent OVS II Opaque and the OVS II Opaque thinner.**

- 18 Dentine opaques, OD10–D41
- 9 Neck opaques
  - OH 10 = Biodent 10
  - OH 13 = Biodent 13, 15, 17, 20, 22
  - OH 16 = Biodent 16
  - OH 21 = Biodent 21
  - OH 23 = Biodent 23, 32
  - OH 25 = Biodent 25, 26, 27, 30
  - OH 31 = Biodent 31

OH 39 = Biodent 39, 41

OH 40 = Biodent 40

4 Intensive opaques

O1 = white

O2 = transparent-grey

O3 = blue-grey

O4 = pink

1 OVS II Opaque thinner

1 set of 3 red sable brushes with metal handles

### Single Components

OVS II Opaque	10 ml
OVS II Opaque thinner	10 ml, 100 ml

## Opaque

The opaque is in liquid form and therefore ready for immediate use. Shake the bottle well before using and replace the stopper immediately after use.

Clean, degrease and dry the metal surface and apply the opaque from the phial with a fine sable brush. The metal must be rough, not smoothed by sand-blasting or electrolytically.

Brush the opaque on thinly and evenly. Apply dentine opaque first, let it dry and then paint on the neck opaque.

With plastic jacket crowns, first insulate the die and then apply opaque to the well-dried insulating film.

## Flask Method

Apply opaque after the crown or bridge has been bedded in and the wax boiled out.

### Drying Time

Opaque hardens at room temperature in 15 minutes. Alternatively, the bridge frame with its covering of opaque may be heated at 60–90°C (e.g. before an infrared lamp) for at least 5 minutes.

## Thin Veneers

The less room there is for the veneer, the harder it is to get the required Biodent shade, depth effect and body. Here the proper use of the additional opaques on top of the neck opaques will be of assistance. Thus, for example, the dentine opaque can be lightened with O1. In the incisal area additional effects can be obtained with opaques O2 and O3.

### Examples of the Use of Opaques O1–O4:

O1 – white: brightens all dentine-coloured opaques

O2 – transparent-grey/O3 – blue-grey:  
To provide more translucency in the incisal area, apply opaque O2 – shaded if necessary with O3 – without perceptibly encroaching on the dentine-coloured opaque, thus avoiding an abrupt transition at the tip.

O4 – pink: covers retentions in contact with the gum

### OVS II Opaque Thinner

If the opaque thickens after a time it can be made less viscous by adding a few drops of thinner (shake well!).

Clean brushes well with the thinner after use – but never use thinner in which brushes have been cleaned for thinning opaques (risk of contamination).

### *Note*

A reliable bond between plastic veneer and metal framework is achieved by the OVS System (Opaque Resin Bonding System). The system consists of coating the metal surface electrolytically with a tin oxide layer to prepare the surface for a micromechanical and chemical bonding with the special OVS II opaque resin (see separate OVS Instruction Manual).

Fig. 3

## **Direct Build-up Method (Liquid S)**

Mix the material to a paste on the palette with the Biodent K+B instrument. The same instrument is suitable for applying and moulding the paste.

Pack undercuts on dies and adjacent plaster teeth with sealing compound and coat the dies and adjacent parts of the cast with Isolation Fluid H. Soak the plaster dies in water before isolating to eliminate air. The isolating film must be left to dry for 15 minutes before applying the plastic paste.

### Mixing

Mix the materials immediately before use. The ratio of liquid and powder is 1:2 by weight.

First pour the monomer into the hollow of the palette and dispense powder until all the monomer is absorbed. The surface is then slightly glossy. Mix well for 15–30 seconds. The proportions are right when the paste hardly drips from the spatula.

Doughing time is from 2 to 4 minutes, depending on room temperature. Then knead the dough well with the spatula. If it is still tacky and cobwebby, moisten the instrument with monomer and smooth the dough. The right consistency is reached when the material is no longer tacky and does not cobweb but is plastic and malleable (if the material is cut into, the cut should not close up). It can now be applied for up to 15 minutes. After that, there is still plenty of time, even in multi-span bridgework, for modelling and carving (another 30 minutes).

### Layering (Figs 4–10)

O = Opaque  
D = Dentine  
S = Enamel  
H = Neck  
T = Transparent

The materials are applied in the order neck, dentine, enamel material.

With thick crowns, an extra wedge of transparent material is laid palatally on the incisal area – with or without an admixture of enamel material (Fig. 4 T). With veneers, the metal frame must first be masked with opaque.

## *Neck material*

Take the work to be veneered, or the insulated die, off the model, if necessary, for layering. Apply the material parallel to the shoulder or to the edge of the gum as shown in the layering diagram (Figs 4–10). Draw it up in the proximal area, as on the shade guide tooth.

In all layers, moisten the instrument with a touch of monomer.

## *Dentine material*

Apply dentine material to the centre of the crown and spread it out. Press it lightly onto the neck area but do not displace the neck material. Undercuts must be packed without forming bubbles. In the incisal area the dentine is built up as with the original Biodent layering (Figs 7–10).

Taper the dentine core off about 1 mm below the incisal edge and bevel off proximally and incisally to a crescent shape. This is the stage at which the dentine core receives its final form.

Fig. 10: Formation of the dentine core.

Blending according to Biodent shade:

Fig. 7: Biodent shades 15, 16, 17

Fig. 8: Biodent shades 21, 23, 31, 39, 41

Fig. 9: Blending for all other Biodent shades.

## *Enamel material*

This is built up in the incisal area and tapered off towards the neck over the dentine.

Take the ready-modelled crown or veneer, with its die, off the model. Carefully remould, or add more material to the edges and points of contact. With veneers, cover the edges with plastic to get a good marginal fit.

## *Protective varnish*

Prevents penetration of the curing medium into the plastic and produces dense, smooth surfaces. **Do not use in cold curing.** Let the work dry after modelling until the surface looks matt. Then brush the varnish on. Leave for about 3 minutes before curing, until the surface is no longer tacky. The drying time may be shortened by warming.

Protective varnish thinner: If the protective varnish thickens after a time, it can be made less viscous again with a few drops of the thinner (shake well).

## Curing

Leave in the vacuum-cure water bath for at least 15 minutes at 95 °C; pressure: 6 bars.

## Fig. 11

Do not exceed the curing temperature of 95 °C.

Cool off slowly, then grind and polish as usual. Fine toothed, hard-metal burs are very suitable for grinding because of the greater hardness.

## Extra Layers

If, perhaps by an error in layering, the dentine shade looks too dark, reduce or roughen the labial surface, moisten with a little monomer and overlay with enamel material. The new material used for correction should have the same consistency as that applied before. Colouring features can now be added or corrections made at the points of contact. The transition from plastic to metal should be masked with fresh material.

## Flask Method (Liquid N)

For veneers, ensure that the walls are thick enough (about 1 mm).

The wax mould should be as close in shape and modelling to the final work as possible, so that extensive grinding may not produce shade variations.

Mechanical retentions are essential in veneering. They should be placed whenever possible in the interdental spaces.

Invest the prepared wax pattern so that the incisal edge rises at an angle of about 30°. The whole labial surface and most of the proximal surfaces should be free of the investment plaster. After boiling-out, coat the plaster mould on both sides with DENTSPLY Isolant.

### Mixing

Mixing ratio of liquid to powder is 1:2 by weight. Pour liquid into the mixing vessel and add enough powder to absorb all the monomer. Mix briefly. If the paste is too thin, add more powder until the paste no longer drips off the spatula. Now mix thoroughly for 15–30 seconds.

A rather higher proportion of liquid is used for the enamel material.

The materials are best mixed in the order in which they are to be used. Doughing time for dentine and neck materials is about 6 minutes. The enamel paste, being thinner, can be applied about 2 1/2 minutes after mixing.

A separate mixing vessel is supplied for each material. The matt surface of the lids can be written on for identification purposes.

### Layering

The materials are layered in the following order: dentine, neck (pre-cure for 10 minutes), enamel (see Figs 4–6). For thick crowns, lay an extra wedge of transparent material – with or without an ad-

mixture of enamel material – palatally on the incisal area (Fig. 4 T). With veneers, the metal frame must first be masked with opaque (see section “Opaque”).

### *Dentine and neck materials*

First mix the dentine material and allow about 6 minutes for doughing. Then pack the hollow mould completely with the paste. Cover with separating foil and close the flask.

Leave the flask in the press and mix the neck material.

Open the flask, cut away the dentine in the neck area and apply the prepared neck material. After a trial closure – and any necessary corrections in the neck area – clamp the flask down and pre-cure for 10 minutes.

### *Enamel material*

After the pre-curing, grind away the dentine in the incisal area as in the Biodent blending diagram (Figs 7–10). Then mix the enamel material, apply it, cover with separating foil and press. Check layering and correct as necessary.

### Curing

Heat the clamped-down flask up to 100 °C in not more than 30 minutes; then boil for 30 minutes more. Cool the flask off slowly to room temperature.

## Cold Curing (Liquid K)

### Applications

- Corrections to crowns and bridge veneers
- Repairs
- Attachment of plastic teeth in the flask-pack method
- Temporary crowns and bridges in the laboratory

Pack undercuts with sealing compound and coat dies and adjacent parts of the cast with Isolation Fluid K. Soak plaster dies in water before isolating to eliminate air. The isolating film must be left to dry for 10 minutes before applying the plastic.

## Mixing

In cold curing the Promotor K is required for hardening. First mix powder and Promotor K in the proportion of 10:1 by weight.

Place 1 large measure of powder and 1 small measure of Promotor K in the special mixing vessel and shake well. Mix this blend on the palette with monomer K to a creamy-to-pasty consistency, depending on the nature of the work.

Working time may be up to 4 minutes after mixing, depending on consistency and room temperature.

The powder/Promotor mixture is workable at room temperature for 4 weeks.

## Layering

As in the direct build-up method described on pages 3–4. Do not use protective varnish!

## Curing

Hardening takes 8–12 minutes from start of mix, depending on room temperature. It is advisable to cure in a pressure vessel for 10 minutes at 40–50 °C.

For large or complicated workpieces an intermediate cure of the applied material is recommended, also of 10 minutes at 40–50 °C, in the pressure vessel.

## Shade Combination Chart

Biodent Shades	Dentine Opaque	Neck Opaque	Dentine Powder	Neck Powder	Enamel Powder
10	OD 10	OH 10	D 10	H 10	S 10
13	OD 13	OH 13	D 13	H 13	S 13
15	OD 15	OH 13	D 15	H 15	S 15
16	OD 16	OH 16	D 16	H 16	S 16
17	OD 17	OH 13	D 17	H 17	S 17
20	OD 20	OH 13	D 20	H 20	S 20
21	OD 21	OH 21	D 21	H 21	S 13
22	OD 22	OH 13	D 22	H 22	S 22
23	OD 23	OH 23	D 23	H 23	S 15
25	OD 25	OH 25	D 25	H 25	S 22
26	OD 26	OH 25	D 26	H 26	S 26
27	OD 27	OH 25	D 27	H 27	S 27
30	OD 30	OH 25	D 30	H 30	S 30
31	OD 31	OH 31	D 31	H 31	S 31
32	OD 32	OH 23	D 32	H 32	S 32
39	OD 39	OH 39	D 39	H 39	S 39
40	OD 40	OH 40	D 40	H 40	S 17
41	OD 41	OH 39	D 41	H 41	S 41

### Nature of Special Risks and Safety Advice

**OVS II Opaque, OVS II Opaque Thinner**  
**Liquids S, N, K**  
**Protective Varnish, Protective Varnish Thinner**

Contains methyl methacrylate.

Do not store above 30 °C.  
 Avoid direct sunlight.

For dental use only.

R37/38

Irritating to respiratory system and skin.

R43

May cause sensitization by skin contact.

S2

Keep out of the reach of children.

S24

Avoid contact with the skin.

S37

Wear suitable gloves.

S46

If swallowed, seek medical advice immediately and show this container or label.

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**DeguDent GmbH**
**DENTSPLY Limited**
**DENTSPLY France SAS**
**DENTSPLY Canada Ltd.**
**DENTSPLY PTY Ltd.**