# D-Light<sup>®</sup> Pro

Instructions for use [EN](#)

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Instrucciones de uso [ES](#)

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Vejledning [DA](#)

Bruksanvisning [SV](#)

Bruksanvisningen [NO](#)

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### **DUAL WAVELENGTH LED CURING LIGHT**

#### **CONTENTS**

1. Safety instructions
2. Features
3. Indications
4. Identification of parts
5. Setup
6. Operation
7. Autoclaving and care of the unit
8. Troubleshooting
9. Disposal
10. Warranty & Repair
11. Packaging
12. Equipment class
13. Specifications
14. Explanation of Symbols and Storage
15. Contact Information
16. EMC Tables

#### **1. SAFETY INSTRUCTIONS**

- The unit should only be used by a qualified dental professional.
- Be sure to follow the instructions in order to use the unit correctly and to safeguard the patient and user.
- Patients as well as dental professionals should always use appropriate safety precautions. Use the D-Light Pro curing light with great care and carefully comply with all safety precautions, including wearing suitable light filtering safety glasses or goggles for the operator, assistant and patient.
- For working with the unit, the owner must provide the person(s) with the written operating instructions in an understandable format as well as in the relevant language(s). The owner is fully responsible to assure that the D-Light Pro unit is in a safe working condition at all times. If in doubt, do not use the unit and contact the supplier.
- Prior to use, make sure that the unit is working normally and safely. It should also be checked on a regular basis.
- Keep the unit out of the reach of infants or children.
- Do not use the unit on patients with a history of eye disease such as cataract or retinal problems. This may cause eye damage.
- The unit must not be disassembled or remodeled. This may cause the unit to leak, generate heat, ignite or explode.
- Only authorized technicians can repair the D-Light Pro curing light, battery packs and charging station.
- Use of accessories not mentioned in this manual may result in unsafe conditions and/or reduced product performance. For this reason, use only factory authorized accessories.
- The unit should not be handled with wet (gloved) hands. This may cause a breakdown.
- The unit should not be exposed to direct sunlight, dust, a wet environment or placed near a heater.
- Make certain there is a match between incoming power supply and the product electrical requirement. Do not use the recharger with voltages different than stated in this manual or labeled on the equipment.
- Do not touch the inside of the battery compartment or the contact points of the battery charger directly with the hands or any metal objects. This may cause a breakdown.
- The battery should not be thrown into a fire or heated. It should not be struck with a sharp object or with force, disassembled or remodeled. This may cause the battery to leak, generate heat, ignite or explode.
- Do not allow the battery contact points to touch other metal objects. Do not carry or store the battery next to metal objects in order to avoid accidental contact with the battery contact points.
- Do not charge, use or leave the battery in high temperatures, near fire or in direct sunlight, etc.
- Store the unit away from flammable materials.
- If the lithium-polymer battery should leak, do not touch the liquid. In case of skin contact, flush immediately with tap water because the liquid may be corrosive and cause skin irritation or damage.
- If any liquid leaks from the battery and comes into contact with the eyes, do not rub. Flush immediately with tap water and seek medical attention. The liquid may cause eye irritation or damage.
- Voltage is present in the charging station and battery packs: do not open! Operate only in dry conditions. Never manipulate with wet hands as this may cause an electrical shock.
- If the battery or unit should begin to smoke, emit a bad odour, deform or discolor during irradiation, charging or storage, then discontinue usage immediately. Continued usage may cause the battery to leak, and/or the unit to generate heat, ignite or explode.
- Avoid sudden or strong impacts to the handpiece. This may cause a breakdown or a reduction in the amount of light emitted.
- Do not use mobile phones or other hand-held electronic devices at the same time as the D-Light Pro.
- The light guide and eye-protection shield (either hard or soft version) should be attached correctly in order to ensure that they do not become detached during use. Check routinely. Should they become loose or cracked, then discontinue their use to avoid hazards such as swallowing or inhalation. Make sure to fix the eye-protection shield correctly.
- If the light guide should crack or break, immediately discontinue use. The light guide is made of glass and there is a possibility of fragments accidentally being displaced into the mouth.
- Do not look directly into the light emitted through the light guide. This may cause eye damage. Prior to use, be sure to fix the eye-protection shield (either hard or soft version).
- Prolonged usage of any curing light can increase temperature. Extended usage near pulp and soft tissues should be avoided to prevent injury. Curing near or around sensitive tissues should be performed in short time

intervals. The Low Power mode (LP) can also be used to limit heat generation when close to the pulp.

- If the handpiece is not used for an extended period of time, or the unit is being transported, then the battery should be removed from the handpiece or charging station.
- Maintain adequate infection control measures when using the D-Light Pro in the oral cavity, and execute the appropriate hygiene plan after usage with patients.
- Protect the D-Light Pro from contamination by using protective plastic barrier-sleeves. These sleeves are intended for single-use; discard after each patient to prevent possible cross-contamination. Protective plastic barrier-sleeves are not needed when the light-guide and handpiece sleeve are autoclaved between each patient.
- When using protective plastic barriers (bags, shields, etc.) with the D-Light Pro, be certain that these are securely affix to avoid patient aspiration and treatment complications.
- Do not use the Detection mode (DT) in conjunction with caries detection dyes.
- Do not use if you have red-green color blindness or visual impairment.

## 2. FEATURES

- LED powered visible light curing and visualisation unit with high power output
- The D-Light Pro can polymerize light cured dental materials, quickly and efficiently.
- The D-Light Pro contains two different types of LEDs with different wavelength peaks. This enables the D-Light Pro to activate commonly used photo initiators used in dentistry, like Camphorquinone (wavelength peak of 468 nm) used in most light-cured materials, and other initiators (wavelength peak of 400 +1-20 nm) used in some light-cured dental materials.
- The D-Light Pro Detection mode (DT) assists the visualisation of bacteria in plaque and carious lesions and the identification of fluorescent materials.
- The D-Light Pro offers an easy-to-use and elegant stainless steel design. The design ensures a user-friendly operation of the unit.
- Lithium polymer batteries, when fully charged, allow for more than 45 individual irradiations of 20 seconds duration each. The lithium polymer battery has no memory effect, and a fast recharging speed of 75 minutes.
- The D-Light Pro is designed to offer a stable and continuous high light output over several years.
- Built-in radiometer
- In addition to the light guide, the handpiece is autoclavable once the internal electronic module and battery pack have been removed.
- The D-Light Pro is under warranty for a period of 3 years, except the battery (1 year warranty) and consumables (no warranty).

## 3. INDICATIONS

- Polymerisation of visible light curing materials curing with a wavelength range of 400 to 480nm. In case of questions about the wavelength ranges for certain materials, please contact their respective manufacturers.
- The Detection mode (DT) assists the visualization of bacteria in plaque, fissures, infected dentin and the presence of micro-leakage. It also helps to identify fluorescent restorative materials, and to evaluate the depth of cracks in the tooth structure.

## 4. IDENTIFICATION OF PARTS

- |                                    |          |
|------------------------------------|----------|
| • Handpiece                        | Figure 1 |
| • Electronic module                | Figure 2 |
| • Battery packs (x2)               | Figure 3 |
| • Light-guide                      | Figure 4 |
| • Soft eye-protector shields (x3)  | Figure 5 |
| • Hard, oval eye-protection shield | Figure 5 |
| • Charging station                 | Figure 6 |
| • Power supply and cable           | Figure 6 |
| • EU/UK plug adapters              | Figure 7 |

## 5. SETUP

- Carefully unpack all D-Light Pro components and familiarize yourself with the individual items. Check for completeness.
- First, check that the voltage indicated on the D-Light Pro charging station name plate complies with the local power supply.

### Battery

- Prior to initial use or use after a long interval of non-use, be sure to charge or recharge the batteries completely. Use the original charging station and lithium polymer battery only. Do not use any other battery chargers or batteries.
- Connect the power cord to the charging station, and the plug of the power cord to the mains power outlet (110-240V AC) (Figure 8)
- Insert the batteries fully onto the docking ports of the charging station until you hear/feel a 'click' (Figure 9)
- There is a LED corresponding to each docking port. When the LED is illuminated, the battery is in the process of being charged. When the LED turns off, then charging is completed. Two batteries can be simultaneously charged.
- Do not connect the depleted battery pack to the charging station until it is properly cleaned and disinfected.
- The time to fully charge an empty battery is approximately 75 minutes.
- When the D-Light Pro is not in use for a prolonged period of time, the battery pack should be disconnected from the handpiece; or the handpiece/battery should be in the "off" position.
- Batteries may be stored on the docking ports of the charging station as long as there is power to the charging station.

### Handpiece

- First, insert the electronic module into the handpiece (Figure 10).
- Next, insert the battery into the rear end of the handpiece/electronic module assembly (Figure 11).

### Light-guide

- Before each use, autoclave the light guide and handpiece.  
**CAUTION!** Be certain to remove the electronic module and battery pack from the handpiece.
- Disinfect the eye-protection shield.
- Insert the light guide into the handpiece (Figure 12). Make sure that the light guide snaps into position. Mount either the hard or soft eye-protection shield on the light guide.
- NOTE: For the Detection mode (DT), the use of the hard eye-protection shield is recommended for improved visibility.

## 6. OPERATION

### Activating the unit

- Gently insert one battery pack into the rear-end of the D-Light Pro handpiece after the electronic module is inserted. You will feel the connector guide the battery automatically into the "Off position".
- When resistance is felt, gently push the battery all the way into the handle. A click will be heard when the battery pack is completely inserted and correctly seated in the D-Light Pro handpiece (Figure 11).
- Rotate the cap of battery pack by quarter turns in either direction, to select the required curing program (Figure 13).
- Adjust the angle of the light if necessary by "rolling" your pen-like grip.
- Press the start button (Figure 14). Keep the tip of the light guide as close as possible to the surface of the material to be cured. Prevent direct contact with the material. In case of contact, carefully remove remaining material from the light guide using a plastic spatula.
- The D-Light Pro handpiece may feel slightly warm during prolonged operation. This is normal.

### Curing programs

| PROGRAM          | FEATURES  |
|------------------|---|
| HP<br>High Power | Standard curing with dual wavelength featuring a 20 second cycle with high output power (approx. 1400 mW/cm <sup>2</sup> ) for maximum efficiency.  |
| LP<br>Low Power  | Low power curing with dual wavelength featuring a 20 second cycle with approx. 700 mW/cm <sup>2</sup> output, optimal for instance when curing close to the pulp.   |
| DT<br>Detection  | Violet light only featuring a 60 second cycle, for the visualization of fluorescent restorative materials and/or identification of bacteria in plaque, fissures, infected dentin, etc. (see Indications). |

### Tone signals

| STONE SIGNAL             | EXPLANATION   |
|--------------------------|---|
| 1x                       | A curing cycle has started or ended.  |
| 1 Quick Tone and a Flash | 10 seconds of curing have been completed.   |
| 2x                       | A curing cycle has been interrupted prematurely.  |
| 3x                       | The internal temperature of the handpiece is too high. After 3 minutes pause, the curing light can be used again. |
| 4x                       | Too many consecutive curing cycles have been performed (>10), and a short pause is required (maximum 3 minutes)   |
| 5x                       | Battery charge is low. The battery needs recharging.  |
| Quick Tone(s)            | Selection of a new program (Battery cap rotation)<br>HP=1 quick tone; LP= 2 quick tones; DT= 3 quick tones        |

### Radiometer Usage

- Select the High Power (HP) curing program
- To confirm proper light output, cover the window of the built-in radiometer in the charging station with the light-guide and press the start button (Figure 15).
- The green LED response means the light output is sufficient for usage.
- The red LED response means that the light output is not sufficient for curing. Check troubleshooting section before contacting technical support.

### 7. AUTOCLAVING AND CARE OF THE UNIT

- The light guide and handpiece sleeve can be steam-autoclaved at 134°C (275°F) maximum (Figure 16).  
**CAUTION!** Prior to autoclaving, remove the electronic module and battery pack from the handpiece by pulling on the "Mode" ring.
- Clean and sterilize the light guide and handpiece in a sterilization bag before each patient application.
- Use a validated sterilization process at a maximum temperature of 134°C and for a time of up to 20 min. Perform sterilization according to EN 17665-1:2006 and EN

556-1:2001 at 134 °C for at least 3 minutes and use steam sterilizers that comply with the requirements of EN 13060:2004+A2:2010, class B or S.

- The eye-protection shields (soft and hard versions) cannot be autoclaved, but should be disinfected using an alcohol-based disinfectant.
- Remove the battery pack prior to routine cleaning, disinfection and maintenance of the unit.
- The docking ports of the charging station can be cleaned with a dry clean brush.
- The charging station, battery, and eye-protection shield should be cleaned with a soft cloth moistened with alcohol. Organic detergents such as thinners or petroleum benzine should not be used. Take care not to get water inside the charging station or on the contact points of the battery.

### 8. TROUBLESHOOTING

If any difficulty is experienced while operating the unit, please check below for possible causes of the problem and suggested corrective actions before seeking further advice or repair from your point of purchase.

| Problem   | Check   | Corrective Action   |
|---|---|---|
| When the start button is pressed, the light will not come on. | Is the battery discharged or the charging status extremely low?                   | Charge the battery.   |
|   | Is the battery set correctly in the handpiece/electronic module assembly?         | If not, reset the battery in the right position (refer to section "activating the unit").   |
|   | Is the electronic module correctly inserted into the handpiece?                   | Remove the electronic module and insert it again into the handpiece   |
|   | Has an extended period of continuous irradiation been carried out?                | The unit is cooling down. After a minimum a 3 minutes where the device is not in use, it can be started again.                                  |
| The battery cannot be charged.                                | Is the charging station correctly connected?                                      | Check to ensure the charging station is correctly connected to the power supply. Check if the power supply is correctly connected to the mains. |
|   | Is the battery firmly set in the charging station?                                | Remove any dust present in the docking ports of the charging station and insert the battery firmly in the docking port.                         |
|   | Is the battery already completely charged?  | Insert the battery into the handpiece and activate to confirm the battery has a charge.   |
| The unit is working normally but the material will not cure.  | Is the shelf life of the material expired?  | Use fresh material.   |
|   | Has the light guide been damaged?   | Clean or replace the light guide.   |
| The light guide or electronic module is difficult to insert   | Is the groove of the light guide or electronic module clean and free from damage? | Apply a small amount of lubricant on the groove of the light guide or electronic module   |
| The radiometer indicates a red light response                 | Has the light guide been damaged?   | Clean or replace the light guide.   |
|   | Is the window of the radiometer free of any material?                             | Clean the radiometer window.  |

## 9. DISPOSAL

Dispose of the battery and all components of the curing light according to local regulations.

## 10. WARRANTY / REPAIR

Warranty: 3 years from the date of purchase for the handpiece, electronic module, charging station and power supply.

In case of a breakdown during the term, repair will be carried out free of charge provided that the unit has been used under normal conditions and according to the instructions for use.

Consumables (such as light guide and eye-protection shields) are not warranted.

The battery is a consumable, but has a 12 month warranty applicable only to battery failure.

In order to benefit from the warranty service, the customer must return the apparatus to be repaired to the GC Europe dealer/importer from which it was purchased, at his own expense.

Before returning the device, please ensure that the products are fully decontaminated and free of debris and other organic materials such as blood residues, saliva and other bodily fluids. Contaminated devices will not be repaired/replaced, and the cost for decontamination of the environment affected by the device will be charged on the basis of the actual cost of the materials and labor required to decontaminate.

The apparatus should be returned suitably packed (possibly in its original packing material), accompanied by all the accessories and by the following information:

- a) Owner's details, including his telephone number.
- b) Details of the dealer/importer.
- c) Photocopy of the consignment note/purchase invoice of the apparatus issued to the owner and indicating, in addition to the date, also the name of the apparatus and its serial number.
- d) A description of the problem.

Transport and any damages caused during transport are not covered by the warranty.

In the event of failures due to accidents or improper use, or if the warranty has lapsed, repairs to the products will be charged on the basis of the actual cost of the materials and labour required for such repairs.

## 11. PACKAGING

### Contents of the kit:

|                                  |   |
|----------------------------------|---|
| Handpiece sleeve                 | 1 |
| Electronic module                | 1 |
| Battery packs                    | 2 |
| 8mm black light-guide (120°)     | 1 |
| Charging station                 | 1 |
| Power supply                     | 1 |
| EU plug adapter                  | 1 |
| UK plug adapter                  | 1 |
| Soft eye-protection shields      | 3 |
| Hard, oval eye-protection shield | 1 |

### Accessories available separately:

Handpiece sleeve, light-guide, battery pack, charging station, power supply with adapters, hard eye-protection shield, soft eye-protection shields (x10)

## 12. EQUIPMENT CLASS

The product complies with all the provisions of the 93/42/EEC Directive on medical devices (as amended by the 2007/47/EC Directive) and of the 2011/65/EU Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment that apply to it.

### Classification of product:

Class I medical device according to Rules 5 and 12 of Annex IX of the 93/42/EEC Directive.

Applied part Type BF



Protection from liquids IP XO

|                   |  |
|-------------------|--|
| EN 980:08         | Symbols for use in the labeling of medical devices   |
| EN 1041:08        | Information supplied by the manufacturer of medical devices  |
| EN 1639:09        | Dentistry - Medical devices for dentistry - Instruments  |
| EN ISO 10650-1:05 | Dentistry - Powered polymerization activators - Part 1: Quartz tungsten halogen lamps  |
| EN ISO 10650-2:07 | Dentistry - Powered polymerization activators - Part 2: Light-emitting diode (LED) lamps   |
| EN ISO 10993-1:09 | Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process   |
| EN ISO 17664:04   | Sterilization of medical devices - Information to be provided by the manufacturer for the processing of resterilizable medical devices   |
| EN 60601-1:05     | Medical electrical equipment - Part 1: General requirements for basic safety and essential performance   |
| EN 60601-1-2:07   | Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance. Collateral standard. Electromagnetic compatibility. Requirements and tests   |
| IEC 60601-2-57:11 | Medical electrical equipment - Part 2-57: Particular requirements for the basic safety and essential performance of non-laser light source equipment intended for therapeutic, diagnostic, monitoring and cosmetic/aesthetic use |
| EN 62471:08       | Photobiological safety of lamps and lamp systems   |

### 13. SPECIFICATIONS

| <b>Handpiece</b>        |   |
|-------------------------|---|
| Light source            | High Power Light Emitting Diode                                   |
| Wavelength              | 400 - 480nm with peak at 400-405nm and 460-465nm                  |
| Average light intensity | 1400 mW/cm <sup>2</sup>   |
| Operation               | Maximum consecutive use, 10 cycles @ 20 sec / 3 minutes pause     |
| Light guide             | Diameter 8mm optical fiber (autoclavable @ 134°C)                 |
| Battery                 | Lithium-polymer, rechargeable, 3.7V, >350mAh                      |
| Battery performance     | >45 cycles @ 20sec  |
| Battery charging        | 75min charging time for a depleted battery pack                   |
| Dimensions              | Diameter: 13-15.2mm<br>Length: 210mm with light guide and battery |
| Weight                  | ~95g  |
| <b>Charging Station</b> |   |
| Input                   | 6 VDC, <1 A   |
| Capacity                | 2 batteries simultaneous charging with overcharging protection    |
| <b>Power Supply</b>     |   |
| Type                    | AC Supply   |
| Input                   | 100-240 VAC, 50/60 Hz, 0.5 A                                      |
| Output                  | 6 VDC / 1 A   |
| Classification          | Class II, Double/reinforced isolated equipment                    |
| Plug                    | Exchangeable wall plug adapters (EU & UK)                         |
| <b>General</b>          |   |
| Environment             | Not intended for use where flammable gases are present            |
| Operating conditions    | 10°C – 30°C   |

#### 14. EXPLANATION OF SYMBOLS AND STORAGE

|   |   |   |  |
|---|---|---|--|
|   | Keep away from sunlight                                   |  | High intensity light   |
|  | Keep dry  |  | Refer to instruction manual / booklet  |
|  | Fragile. Handle with care.                                |  | Attention, see instructions for use  |
|   | Temperature limitations:<br>-20°C to +45°C                |  | Please do not discard this device in household garbage.<br>(See "Disposal" section above)  |
|   | Humidity limitation:<br>10 % to 95 %<br>relative humidity |  | Handpiece sleeve & light-guide:<br>Sterilizable in a steam sterilizer (autoclave) at 134°C |
|   | Atmospheric pressure limitations:<br>50 kPa to 106 kPa    |  | CE marking of conformity   |

#### 15. CONTACT INFORMATION

In case of any questions, please contact the GC Europe dealer/importer from which the product was purchased.

| Manufacturer  |   |
|---|---|
| GC Europe N.V.<br>Researchpark Haasrode-Leuven 1240<br>Interleuvenlaan 33<br>B-3001 Leuven, Belgium |  |

#### 16. EMC TABLES

The D-Light Pro is intended for use in the electromagnetic environment specified below. The customer or the user of the D-Light Pro should ensure that it is used in such an environment.

##### Guidance and manufacturer's declaration – Electromagnetic emissions

EN 60601 – 1-2 / Tabel 201

|   |          |  |
|---|----------|--|
| RF emissions CISPR 11                                     | Group 1  | The D-Light Pro uses RF energy only for its internal function. Therefore its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  |
| RF emissions CISPR 11                                     | Class B  |  |
| Harmonic emissions<br>IEC 61000-3-2                       | Complies | The D-Light Pro is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Voltage fluctuations / flicker<br>emissions IEC 61000-3-3 | Complies |  |

**Guidance and manufacturer's declaration – Electromagnetic immunity**

EN 60601 – 1-2 / Table 202

| Immunity test  | Test level   | Compliance   | Electromagnetic environment guidance   |
|--|--|--|--|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2   | ±6 kV contact<br>±8 kV air                                     | The device continues to work regularly and in safety.  | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.  |
| Radiated RF<br>IEC 61000-4-3   | 3 V/m 80 MHz to 2.5 GHz  | The device continues to work regularly and in safety   | Field strengths from fixed transmitters cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered.  |
| Electrical fast transient / burst<br>IEC 61000-4-4   | ±2 kV for power supply lines<br>±1 kV for input / output lines | The device continues to work regularly and in safety   | Mains power quality should be that of a typical commercial or hospital environment.  |
| Surge<br>IEC 610004-5  | ±1 kV differential mode<br>±2 kV common mode                   | The device continues to work regularly and in safety.  |  |
| Conducted RF<br>IEC 61000-4-6  | 3 Veff 150 kHz to 80 MHz                                       | The device continues to work regularly and in safety.  | Portable and mobile RF communications equipment should be used no closer to any part of the disposal, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br>Recommended separation distance1):<br>$d = 1.2 \sqrt{P}$<br>$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz<br>$d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000   | 3-100 A/m  | The device continues to work regularly and in safety.  | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.  |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 | <5% U/10ms<br>70% U/ 0.5s<br>40% U/0.1s                        | The device can vary from the required levels of immunity with a duration for as long as the device remains in safety; no malfunctions have been detected and can be restored to pre-test status with the intervention of the operator. | Mains power quality should be that of a typical commercial or hospital environment.  |

**Notes:**

- 1.) (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and (d) is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol: 
- 2.) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the D-Light Pro is used exceeds the applicable RF compliance level above, the D-Light Pro should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the D-Light Pro.

**Recommended separation distances between portable and mobile RF communications equipment and the D-Light Pro**

EN 60601-1-2 - Table 206

This product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the D-Light Pro can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the D-Light Pro as recommended below, according to the maximum output power of the communications equipment.

| Separation distance according to frequency of transmitter (m) |   |  |   |
|---|---|--|---|
| Rated maximum output power of transmitter (W)                 | 50 kHz-80 MHz<br>$d = 1,2 \times \sqrt{P}$ Sender | 80 MHz-800 MHz<br>$d = 1,2 \times \sqrt{P}$ Sender | 800 MHz-2,5 GHz<br>$d = 2,3 \times \sqrt{P}$ Sender |
| 0,01  | 0,12  | 0,12   | 0,23  |
| 0,1   | 0,38  | 0,38   | 0,73  |
| 1   | 1,2   | 1,2  | 2,3   |
| 2   | 1,7   | 1,7  | 3,25  |
| 10  | 3,8   | 3,8  | 7,3   |
| 100   | 12  | 12   | 23  |

## IFU ILLUSTRATIONS

| Figures   | Captions   |
|-----------|--|
| Figure 1  | <b>Handpiece</b>   |
| Figure 2  | <b>Electronic module</b>   |
| Figure 3  | <b>Battery packs (x2)</b>  |
| Figure 4  | <b>Light-guide</b>   |
| Figure 5  | <b>Eye-protection shields</b><br>Slide the eye-protection shield (soft or hard version) on the end of the light guide  |
| Figure 6  | <b>Charging station and power supply</b>   |
| Figure 7  | <b>Attaching the plug adapter</b>  |
| Figure 8  | <b>Connecting the charging station</b>   |
| Figure 9  | <b>Recharging the battery packs on the charging station</b>  |
| Figure 10 | <b>Inserting the electronic module into the handpiece</b>  |
| Figure 11 | <b>Inserting the battery pack into the handpiece</b>   |
| Figure 12 | <b>Inserting the light guide into the handpiece</b><br>Adjust the light angle; the light guide rotates 360°<br>Position the light guide tip close to the material for best results<br>Check the light guide attachment regularly<br>Check and clean both ends of the light guide   |
| Figure 13 | <b>Selecting the curing program</b><br>Rotate the battery pack by quarter turns until your desired program faces the reference point.  |
| Figure 14 | <b>Activating the curing light</b><br>Press the switch on the handpiece.<br>The curing light is activated for 20 seconds. After 10 seconds have passed, there will be a quick tone and a flash of the light. At the end of the curing cycle, a tone will sound and the curing light will stop. You can also press the switch to interrupt the curing cycle before its end. |
| Figure 15 | <b>Using the radiometer</b><br>Position the light-guide as shown in the illustration and activate the light in High Power mode (HP). A green light indicates sufficient power output, while a red light warns of an insufficient power.  |
| Figure 16 | <b>Autoclaving the handpiece sleeve &amp; light-guide</b><br><b>CAUTION!</b> Make sure to remove the electronic module and battery pack before placing the handpiece in the autoclave.   |

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3737 West 127th Street, Alsip, IL 60803 U.S.A.
**GC ASIA DENTAL PTE.LTD**
**GC AUSTRALASIA DENTAL PTY.LTD**