
# GC FUJIVEST® PLATINUM II

|    |                                                        |
|----|--------------------------------------------------------|
| EN | PHOSPHATE BONDED INVESTMENT (TYPE I) CLASS 1&2         |
| DE | PHOSPHATGEBUNDENE EINBETTMASSE (TYP I) KLASSE 1 UND 2  |
| FR | REVETEMENT A LIANT PHOSPHATE (TYPE I) CLASSE 1&2       |
| IT | RIVESTIMENTO A LEGANTE FOSFATICO (TIPO I) CLASSE 1 E 2 |
| NL | FOSFAATGEBONDEN INBEDMASSA (TYPE I) KLASSE 1 EN 2      |
| ES | REVESTIMIENTO DE BASE DE FOSFATO (TIPO I) CLASE 1 Y 2  |
| SV | FOSFATBUNDEN INBÄDDNINGSMASSA (TYP I) KLASSE 1&2       |
| DA | FOSFATBUNDEN INDSTØBNINGSMASSE (TYPE I) KLASSE 1&2     |
| FI | FOSFAATTISIDONNAINEN VALUMASSA (TYYPPI I) LUOKKA 1&2   |

|    |                                                         |
|----|---------------------------------------------------------|
| EN | NEW LIQUID DILUTION CHART – SEE PAGE 5                  |
| DE | NEUE TABELLE ZUR FLÜSSIGKEITSVERDÜNNUNG – SIEHE SEITE 8 |
| FR | NOUVEAU TABLEAU DE DILUTION DU LIQUIDE – PAGE 11        |
| IT | NUOVA TABELLA DILUIZIONE – VEDI PAG. 14                 |
| NL | NIEUWE VLOEISTOFVERDUNNINGSKAART – ZIE PAGINA 17        |
| ES | NUEVA TABLA DE DILUCIÓN LÍQUIDO - VER PÁGINA 20         |
| SV | NYTT SCHEMA FÖR VÄTSKEUTSPÄDNING – SE SIDAN 23          |
| DA | NEW LIQUID BLANDINGSFORHOLD – SE SIDE 26                |
| FI | UUSI NESTEEN LAIMENNUSTAULUKKO – KATSO SIVU 29          |

**FEBRUARY 2018**

**GC**

# GC FUJIVEST® PLATINUM II

|    |                                  |    |
|----|----------------------------------|----|
| EN | PHOSPHATE BONDED INVESTMENT      | 4  |
| DE | PHOSPHATGEBUNDENE EINBETTMASSE   | 7  |
| FR | REVETEMENT A LIANT PHOSPHATE     | 10 |
| IT | RIVESTIMENTO A LEGANTE FOSFATICO | 13 |
| NL | FOSFAATGEBONDEN INBEDMASSA       | 16 |
| ES | REVESTIMIENTO DE BASE DE FOSFATO | 19 |
| SV | FOSFATBUNDEN INBÄDDNINGSMASSA    | 22 |
| DA | FOSFATBUNDEN INDSTØBNINGSMASSE   | 25 |
| FI | FOSFAATTISIDONNAINEN VALUMASSA   | 28 |

Dear customer,

We thank you for purchasing a high quality GC product.

This Instruction For Use gives you a clear summary of each working stage and is based on the results of numerous lab tests carried out by our Research and Development Department.

If you have any further questions on the use of this product, please contact your local GC representative.

For more information on GC products please visit our website [www.gcamerica.com](http://www.gcamerica.com).

#### GENERAL REMARK

All information included in this Instruction For Use is based on extensive testing and an extended range of casting trials. However, because of different working methods and equipment (e.g. waxes, resins, casting liner, mixing equipment, etc.), different end-results may be obtained.

GC FUJIVEST® PLATINUM II is a carbon-free phosphate bonded precision casting investment for quick heating as well as conventional step heating procedure.


Specially developed for precision castings of precious, semi-precious and Pd-base alloys for crown and bridgework.

GC FUJIVEST® PLATINUM II offers you a number of advantages:

- Can be used in a quick-heating method or a progressive heating-up schedule
- Extremely smooth casting surfaces
- Perfect fit of the castings even for long span restorations like implant works
- Very easy devesting
- Easy mixing because of optimal wettability

PHYSICAL DATA (typical data)

|                                                                      |                      |                                  |
|----------------------------------------------------------------------|----------------------|----------------------------------|
| ISO 15912: 2006 ; Measured at 23°C / 73°F with 100% Expansion liquid | 20' setting          | 120' setting                     |
| <b>Heating up / Burn out procedure</b>                               | <b>Quick heating</b> | <b>Conventional step heating</b> |
| Setting expansion                                                    | 0.70%                | 1.70%                            |
| Thermal expansion                                                    | 1.40%                | 0.80%                            |
| Total expansion (Linear)                                             | 2.10%                | 2.50%                            |
| Initial setting time                                                 | 7 minutes            | 7 minutes                        |
| Compressive strength                                                 | 5 Mpa                | 5 Mpa                            |
| Working time                                                         | 6 minutes            | 6 minutes                        |
| Flow                                                                 | 62,5 mm              | 62,5 mm                          |

## INSTRUCTIONS FOR USE

### 1. Storage

Recommended for optimal performance store powder and liquid at normal room temperature (23°C / 73°F).

In case the powder was exposed to lower temperature, it is necessary to equilibrate the material at room temperature for at least 24 hours before use. Storage and working temperature of investment powder and liquid are an important factor in determining the setting time and expansion, and hence the fit and the surface roughness of the castings.

- ① Protect the liquid against low temperatures!  
Do not store below 23°C / 73°F.  
Once frozen, the liquid cannot be used anymore.  
Be careful with winter deliveries.

### 2. Working temperature

Use at ± 23°C / 73°F room temperature (19°C minimum).

- ① Higher working temperatures slightly reduce working time but improve the surface smoothness.  
Lower temperatures prolong setting time and could result in surface roughness.

### 3. Preparations before investing

During wax modeling use GC Multi Sep as a wax separator, which leaves no oily film on the die surface.

Use a 1 mm thick dry casting liner, Casting Ring Liner.

Seal the borders of the liner with a thin layer of vaseline.

- ① Open or wet liners can change the liquid content of the mixture and thus affect the expansion.  
Be sure that the total inner surface of the metal ring is covered with casting liner.  
Ensure that the liner is flush with the upper rim of the metal ring, so no direct contact of investment with the metal ring.  
X3 ringsize = 1 layer of Casting Ring Liner  
X6 ringsize = 2 layers of Casting Ring Liner  
X9 ringsize = 2 layers of Casting Ring Liner

GC FUJIVEST® PLATINUM II is a very fluid investment, which can be used without any wetting agents.

Nevertheless, Wax Modellations may be appropriately treated with a surface tension reducing liquid/spray in order to permit better adherence and flow of investment.

If a surface agent is used, be sure to totally dry the surface prior to investing.

Too much surface agent or choice of inappropriate surface agent could cause rough surface or investment residues in the casted object.

### 4. Ring sizes

Can be used with a ring and ringless technique.

- ① Possible to use with: X1, X3, X6, X9 for ring technique and X1, X3, X6 for ringless technique
- ① With ringsize X1 care should be taken. Due to smaller amount of material, meaning less exothermic reaction during setting, smaller castings could result. Temperatures below 23°C / 73°F may result in rough surfaces. Ringsize X3 gives the best results.

## 5. Powder/liquid ratio

| Ringsize | Powder (g) | Liquid (mL) |
|----------|------------|-------------|
| x1       | 60         | 13,2        |
| x3       | 150        | 33          |
| x6       | 300        | 66          |
| x9       | 450        | 99          |

- ① Standard P/L ratio = 100 g / 22 mL.  
Exact powder/liquid measurement is necessary to obtain stable results.  
Use adapted measure equipment such as electric balance, liquid measuring cup, pipette.

## 6. Expansion

Basic info on liquid concentration.

The used concentration of the expansion liquid depends on various parameters:

- storage temperature of the investment (powder and liquid)
- mixing time
- setting under pressure
- the amount of mixture

- ① Pure liquid gives the maximum setting expansion, meaning the largest casting.  
Pure water gives the minimum setting expansion, meaning the smallest casting.  
Use only distilled water to dilute.

## Liquid dilution chart

| Type Alloy                          |                        | Ring-size | High Precious<br>>70% Au |                | Semi Precious<br>< 55% Au |                | Pd -Base  |                | High Precious<br>ceramic alloy |                | Semi Precious<br>ceramic alloy |                | Pd-Base<br>ceramic alloy |                |
|-------------------------------------|------------------------|-----------|--------------------------|----------------|---------------------------|----------------|-----------|----------------|--------------------------------|----------------|--------------------------------|----------------|--------------------------|----------------|
|                                     |                        |           | mL liquid                | mL dist. water | mL liquid                 | mL dist. water | mL liquid | mL dist. water | mL liquid                      | mL dist. water | mL liquid                      | mL dist. water | mL liquid                | mL dist. water |
| Inlays / Onlays<br>Partial Crowns   | Wax Pattern            | x1        | 4,4                      | 8,8            | 4,4                       | 8,8            | 6,8       | 6,4            |                                |                |                                |                |                          |                |
|                                     |                        | x3        | 11                       | 22             | 11                        | 22             | 17        | 16             |                                |                |                                |                |                          |                |
|                                     |                        | x6        | 22                       | 44             | 22                        | 44             | 34        | 32             |                                |                |                                |                |                          |                |
|                                     |                        | x9        | 33                       | 66             | 33                        | 66             | 51        | 48             |                                |                |                                |                |                          |                |
| Crowns<br>Bridges<br>Primary Crowns | Wax Pattern            | x1        | 4,8                      | 8,4            | 5,6                       | 7,6            | 7,6       | 5,6            | 7,6                            | 5,6            | 6,8                            | 6,4            | 6,8                      | 6,4            |
|                                     |                        | x3        | 12                       | 21             | 14                        | 19             | 19        | 14             | 19                             | 14             | 17                             | 16             | 17                       | 16             |
|                                     |                        | x6        | 24                       | 42             | 28                        | 38             | 38        | 28             | 38                             | 28             | 34                             | 32             | 34                       | 32             |
|                                     |                        | x9        | 36                       | 63             | 42                        | 57             | 57        | 42             | 57                             | 42             | 51                             | 48             | 51                       | 48             |
|                                     | GC Pattern<br>Resin LS | x1        | 6                        | 7,2            | 6,4                       | 6,8            | 8         | 5,2            |                                |                |                                |                |                          |                |
|                                     |                        | x3        | 15                       | 18             | 16                        | 17             | 20        | 13             |                                |                |                                |                |                          |                |
|                                     |                        | x6        | 30                       | 36             | 32                        | 34             | 40        | 26             |                                |                |                                |                |                          |                |
|                                     |                        | x9        | 45                       | 54             | 48                        | 51             | 60        | 39             |                                |                |                                |                |                          |                |
| Secondary<br>Units                  | GC Pattern<br>Resin LS | x1        | 4,8                      | 8,4            | 5,6                       | 7,6            | 6,8       | 6,4            | 6,8                            | 6,4            | 6                              | 7,2            | 6                        | 7,2            |
|                                     |                        | x3        | 12                       | 21             | 14                        | 19             | 17        | 16             | 17                             | 16             | 15                             | 18             | 15                       | 18             |
|                                     |                        | x6        | 24                       | 42             | 28                        | 38             | 34        | 32             | 34                             | 32             | 30                             | 36             | 30                       | 36             |
|                                     |                        | x9        | 36                       | 63             | 42                        | 57             | 51        | 48             | 51                             | 48             | 45                             | 54             | 45                       | 54             |

## NOTE

All information included in this Instruction For Use is based on extensive testing and a whole series of casting trials.

However, because of different working methods and equipment (e.g. waxes, resins, casting liner, mixing equipment, etc.) different end-results may be obtained.

All fit tests have been carried out with the quick heating method (20 minute setting)

## 7. Mixing

1. Pre-mix powder and liquid thoroughly by hand with a spatula.

- ① Be sure all powder is wetted out by the liquid to give a uniform mixture.

2. Mix for 60 seconds under vacuum (420 rpm).

- ① Always use clean mixing bowl and check vacuum level. Insufficient vacuum leads to differences in fit and air-bubbles in the castings.

## 8. Working time

At least 4 minutes pouring time at room temperature (23°C / 73°F).

- ① Working time / pouring time depends on the temperature of powder & liquid and on room temperature. Higher temperatures shorten the working time.

## 9. Investing

Invest under gentle vibration (low frequency).

- ① GC FUJIVEST® PLATINUM II is very fluid, so strong vibration is not necessary and not advisable.  
From the moment the ring is totally filled, **stop vibration immediately** and do not touch the investment until set.

## 10. Setting time

Leave to set for 20 minutes from **start of mixing**.

- ① **Best results are obtained by putting immediately into a preheated furnace after 20 minutes.**

Scrape the top surface of the investment ring with a sharp knife.

- ① Diamond disks for trimmers are damaged by the quartz and cristobalite particles of the investment.

## 11. Heating Up – Burn Out Procedure

|                       | Quick heating                                                                                                  | Conventional Step heating                                                                                                                                                                          |
|-----------------------|----------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Insertion temperature | Pre-heated furnace at<br>750°C-800°C/1380°F-1470°F for Au-alloys<br>800-850°C*/1470-1560°F* for Ceramic alloys | Room temperature                                                                                                                                                                                   |
| Step 1                |                                                                                                                | Room temperature (23°C/73°F) to 260°C/500°F<br>Heat rate 3°C/37°F minutes                                                                                                                          |
| Step 2                |                                                                                                                | Holding time at 260°C/500°F<br>60 - 90 minutes                                                                                                                                                     |
| Step 3                |                                                                                                                | Rise temperature from 260°C/500°F to 580°C/1076°F<br>5°C/41°F per minute                                                                                                                           |
| Step 4                |                                                                                                                | Holding time at 580°C/1076°F<br>20 - 50 minutes                                                                                                                                                    |
| Step 5                |                                                                                                                | Rise temperature from 580°C/1076°F to 750°C-800°C/1380°F-1470°F for Au-alloys<br>Rise temperature from 580°C/1076°F to 800-850°C*/1470-1560°F* for Ceramic alloys<br>Heat rate 5°C/41°F per minute |
| Holding time          | X 1 40 minutes at end temperature                                                                              | X 1 30 minutes at end temperature                                                                                                                                                                  |
|                       | X 3 50 minutes at end temperature                                                                              | X 3 40 minutes at end temperature                                                                                                                                                                  |
|                       | X 6 60 minutes at end temperature                                                                              | X 6 50 minutes at end temperature                                                                                                                                                                  |
|                       | X 9 90 minutes at end temperature                                                                              | X 9 60 minutes at end temperature                                                                                                                                                                  |

Due to aggressive burning out, do not open the furnace during heating-up.

In case of vacuum pressure casting, raise the end temperature by 50°C/122°F.

\* 850°C/1560°F: for casting-on technique of prefabricated metal parts.

When several investment rings are put into the furnace at the same time, prolong the heating period.

**750°C/1380°F is the minimum pre-heating temperature**

## 12. Casting

Cast in the usual manner: centrifugal casting, vacuum-pressure casting, etc.

Take care with the positioning of the casting ring in the casting device.

Cast as soon as possible after removing the ring from the furnace.

## 13. Cooling

Cool down the casting as slowly as possible.

e.g. insert into cold furnace and close the door.

Place the ring upside down for cooling down.

## Notes

1. Clean bowl, spatula, etc. carefully after use. The chemical composition of the residual GC FUJIVEST® PLATINUM II will delay the setting time of gypsum products.
2. It is recommended that mixing bowls are stored in water between uses.
3. Store powder & liquid at room temperature (+/- 23°C / 73°F).

## Safety recommendations and hazard warnings

1. Before use read the powder and liquid products Safety Data Sheets communicated by your dealer.  
These are also available on <http://www.gceurope.com> or for the US and Canada from <http://www.gcamerica.com>
2. The solid product contains quartz and cristobalite. Causes damage to the lung through prolonged or repeated exposure.  
Route of exposure: Inhalative. Do not breathe dust. In case of inadequate ventilation wear respiratory protection.
3. Open the investment material bag with scissors and avoid the formation of dust when filling into the mixing bowl.  
Rinse the empty investment material bag with water before disposal.
4. Remove dust from your working place only when it is wet.
5. To avoid the formation of dust when removing the investment material from the casting ring,  
place the cooled casting ring into water for a short time.
6. When sandblasting the cast object, always use a fine-dust filter extraction system.
7. During heating up the material, use of a fume hood is mandatory. Working in closed places is to be avoided.  
Do not inhale fumes when heating the material.
8. Never touch hot materials by hand, always use appropriate tools like muffle grippers in order to avoid burns.
9. Take care of the heat caused by the chemical setting of the material, in order to avoid harm.
10. Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

## Packages

Powder: 6 kg box (60 g pack x 100)  
6 kg box (90 g pack x 67)  
6 kg box (150 g pack x 40)

Liquid: 900 ml bottle

Expiry date: 2 years from the manufacturing date

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**GC NORDIC AB**