
R-030350

ES..... PRENSA HIDRÁULICA	7
EN..... HYDRAULIC PRESS	8
FR..... PRESSE HYDRAULIQUE	9
PO PRENSA HIDRÁULICA	10
IT..... PRESSA IDRAULICA.....	11
DE..... HYDRAULIKPRESSE	12

MESTRA®

TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tel. +34 944530388 - Fax +34 944711725

48150 SONDIKA - BILBAO - ESPAÑA

mestra@mestra.es - www.mestra.es

Rev. 13/12/17


Fig. 1
Abb.1

- | | |
|---|--|
| <p>A Husillo de aproximación / Approach spindle / <i>Vis axiale</i> / Fuso de aproximação / Vite di avvicinamento / <i>Gewindespindel</i></p> <p>B Tornillo de cierre / Locking screw / <i>Vis de blocage</i> / Parafuso de fecho / Vite di chiusura / <i>Verschlusschraube</i></p> | <p>C Tornillo de carga de aceite / Oil filling screw / <i>Vis-bouchon de remplissage d'huile</i> / Parafuso de carga de óleo / Vite di carico olio / <i>Öfüllschraube</i></p> <p>D Palanca / Lever / <i>Palan</i> / Alavanca / Leva / <i>Hebel</i></p> <p>E Platinos superior e inferior / Top & bottom plates / <i>Plateaux (supérieur et inférieur)</i> / Discos superior e inferior / Piatti superiore e inferiore / <i>Teller und Hubtisch</i></p> |
|---|--|

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This ring is only for transportation. Do not try to relocate it.

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


SAFETY MEASURES FOR THE R-030350 HYDRAULIC PRESS

- NEVER EXCEED THE RATED CAPACITY OF THE PRESS.
- DO NOT PUT YOUR HANDS OR ANY OTHER PART OF YOUR BODY INTO THE OPERATIONS AREA.
- MAKE SURE THE PART IS POSITIONED CORRECTLY BEFORE COMMENCING THE OPERATION.
- FAILURE TO COMPLY WITH THESE SAFETY MEASURES MAY CAUSE DAMAGE OR INJURY TO THE USER, THE PRESS OR THE PART.

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The R-030350 is a practical hydraulic press, with maximum power level 15.2 tonnes (rated, 9.5 tonnes), specially designed for use in the manufacture of muffles for dental prostheses.

INSTRUCTIONS FOR USE

1. Put the muffle(s) between the plates (E) taking care to centre it/them properly.
2. Bring up spindle (A) until it touches the muffle(s), close up the screw (B) and start pumping.
3. The pressure gauge will not stabilize until the muffle(s) is/are compacted and all excess resin is squeezed out. When the pressure gauge needle is steady, optimum pressing has been achieved.
4. To remove the muffle(s) from the press, unscrew the locking screw (B) so that the bottom plate (E) is lowered by the springs, and the unit will be ready for the next pressing.

Every 50 Bars (725 psi) a pressure of 1,900 kg is obtained. Thus at 100 bars (1450 psi) 3,800 kg is obtained, and at 200 bars (2900 psi) 7,600 kg is obtained. The maximum is 400 bars (5801 psi) or 15,200 kg of pressure. We do not advise going this far, as the pressure gauge may suffer mechanical damage at this level, and the service life of all components of the press and the muffles may be reduced.

PRECAUTIONS

- ◆ Do not operate the press with fragile materials.
- ◆ Protect yourself adequately during operation of the unit, and use safety gauntlets and goggles.
- ◆ Do not allow children or unqualified staff to handle the unit.
- ◆ Before using the facility, ensure that it is resting on a firm solid base.
- ◆ Ensure that the items to be pressed are correctly positioned between the plates.
- ◆ Be extremely careful to remove your fingers from the plates before starting the pressing operation.

TECHNICAL CHARACTERISTICS

Height:	510 mm
Width:	225 mm
Depth:	230 mm
Table diameter:	125 mm
Effective height:	235 mm
Effective width:	140 mm
Maximum piston stroke:	25 mm
Weight:	22 kg
Maximum power permitted:	400 bars (5801 psi) → 15200 kg
Maximum power recommended:	250 bars (3625 psi) → 9500 kg

The maximum recommended pressure is 250 bars (3625 psi), red mark on pressure gauge. Normal working pressure is 150 - 250 bars (2175 – 3625 psi). If more liquid is needed use only EXTRA (SAE 10) oil. No other liquids may be used under any circumstances.

IMPORTANT

During the warranty period this press should only be handled at our workshops. Any handling elsewhere will cause cancellation of the warranty.

All assembly and disassembly operations should be carried out with extreme care so as not damage any part of the press.

If there is considerable oil spillage between the piston of the pump (3) and the nut (4) when the lower plate drops to its rest position, too much oil has been put into the unit (Fig. 1). Unscrew the oil filling screw and allow excess oil to drain off. The level must be between Max. and Min. (Fig. 1).

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