
# Autopol *Geisser*

Ref. 030425

ES.....	7	Manual original
EN.....	12	
FR.....	27	
DE.....	22	


N.º	Ref.	DESCRIPCIÓN	DESCRIPTION
1	030425-25	Manguera teflón Ø 6x8	Teflon tube Ø 6x8
2	080561-47	Manguera teflón Ø 4x6	Teflon tube Ø 4x6
3	000250	Cable	Cable
4	000231	Fusible 4 A	Fuse 4 A
5	080490-06	Base de conexión	Power socket
6	080422-17	Interruptor	Switch
7	080518-11	Filtro electroválvula	Electro valve filter
8	030425-43	Racor entrada 1/8 Ø6	Inlet fitting 1/8 Ø6
9	030425-33	Racor recto 1/8 Ø8	Straight fitting 1/8 Ø8
10	030425-05	Chapa trasera	Back cover
11	030425-34	Racor en T	T fitting
12	030425-35	Racor alargador	Long fitting
13	030425-36	Válvula salida agua	Water outlet valve
14	030425-19	Racor recto 1/8 Ø6	Straight fitting 1/8 Ø6
15	030425-37	Tapón	Cap
16	080253-10	Electroválvula 24 V	Solenoid valve 24 V
17	080460-69	Racor botella 1/8 Ø8	Bottle fitting 1/8 Ø8
18	080460-89	Tuerca	Nut
19	080460-77	Bidón	Bottle
20	030425-21	Racor HH 1/8	Fitting FF 1/8
21	030425-39	Espiga	Peg
22	030425-40	Regulador presión	Pressure regulator
23	030425-41	Racor MM 1/8	Fitting MM 1/8
24	030425-52	Válvula seguridad 2.2 bar	Safe valve 2.2 bar
25	030425-42	Racor 4 vías	4-way fitting
26	030425-20	Racor codo 1/8 Ø6	Elbow fitting 1/8 Ø6
27	20-080614	Fuente alimentación 24 V	Power source 24 V
28	080508-18	Pata	Leg

N.º	Ref.	DESCRIPCIÓN	DESCRIPTION
29	030425-03	Carcasa inferior	Inner housing
30	030425-06	Soporte olla	Tank support
31	030425-10	Resistencia	Heating element
32	030425-02	Carcasa superior	Upper housing
33	100210-08	Junta carcasa	Housing gasket
34	030425-09	Manómetro	Pressure gauge
35	030425-27	Junta tórica manómetro	Pressure gauge O-ring
36	030425-22	Codo manómetro	Pressure gauge elbow
37	030425-51	Tuerca ciega	Cap nut
38	030425-15	Carátula control	Control cover
39	030425-11	Control electrónico	Electronic control
40	080511-63	Junta tórica sonda	Thermal probe O-ring
41	030425-12	Sonda temperatura	Thermal probe
42	030425-44	Alargador racor distribuidor	Distributor fitting extension
43	030425-14	Grapa tubo aire	Air tube clip
44	030425-01	Cuba	Tank
45	030425-04	Embellecedor tapa	Lid cover
46	030425-52	Válvula seguridad 2.2 bar	Safe valve 2.2 bar
47	030420-13	Arandela cobre	Copper ring
48	030425-26	Tuerca válvula tapa	Lid valve nut
49	030425-30	Junta tapa	Lid ring
50	030425-13	Bandeja	Tray
51	030425-07	Distribuidor central	Central distributor
52	030425-08	Junta tórica 63 x 2	O-ring 63 x 2
53	030425-17	Clixon 140 °C	Klixon 140 °C
54	030425-54	Asa con pulsador	Handle with button
55	030425-55	Asa sin pulsador	Handle without button
56	030425-45	Accesorio: adaptador 1.5 bar	Accessory: 1.5 bar adaptor

F

G

H

I

J

K

L Ref. 030425-45

- A Tapa / **Lid** / *Couvercle* / Deckel
  - B Olla / **Tank** / *Réceptient* / Behälter
  - C Teclado / **Keyboard** / *Clavier* / display mit Bedienungstasten
  - D Apertura tapa / **Lid opening** / *Ouverture du couvercle* / Deckelöffnung
  - E Manómetro / **Pressure gauge** / *Manomètre* / Manometeranzeige
  - F Llave agua / **Water valve** / *Valve pour l'eau* / Kugelhahn
  - G Salida agua / **Water outlet** / *Sortie eau* / Wasserablauf
  - H Entrada aire / **Air inlet** / *Entrée de l'air* / Lufteinlass
  - I Salida aire / **Air outlet** / *Sortie air* / Luftauslass
  - J Conexión eléctrica / **Power supply** / *Connexion électrique* / Stromanschluss
  - K Interruptor general / **Main switch** / *Interrupteur général* / Hauptschalter
- 
- L Accesorio: Adaptador 1.5 bar / **Accessory not included: 1.5 bar adaptor** / *Accessoire non inclus: Adapteur 1.5 bar* / Zubehör nicht enthalten: 1.5 bar

**TALLERES MESTRAITUA, S.L.**  
Txori-erri Etorbidea, 60  
Tel. + 34 944 530 388\* - Fax + 34 944 711 725  
mestra@mestra.es - www.mestra.es  
**48150 SONDIKA - BILBAO - (España)**

**TALLERES MESTRAITUA, S.L.**  
Txori-erri Etorbidea, 60  
Tel. + 34 944 530 388\* - Fax + 34 944 711 725  
mestra@mestra.es - www.mestra.es  
**48150 SONDIKA - BILBAO - (España)**

**TALLERES MESTRAITUA, S.L.**  
Txori-erri Etorbidea, 60  
Tel. + 34 944 530 388\* - Fax + 34 944 711 725  
mestra@mestra.es - www.mestra.es  
**48150 SONDIKA - BILBAO - (España)**

**GARANTIA**

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**GUARANTEE**

All MESTRA® units are guaranteed against every defect in manufacturing for a period of one year from date of purchase. Prior to honouring a guarantee, MESTRA® will require that the attached coupon be filled in with date of sale, duly stamped by the distributor.

The guarantee does not cover breakdowns caused by incorrect installation or use.

MESTRA® will not effect any repair outside its own workshops and will require transport costs.

Repairs carried out on the premises of the customer, be they with or without guarantee, are entrusted to the technical assistance of the distributor who has sold the unit or alternatively to the technicians he has contracted

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

Director Gerente

Sondika, 08.07.2019

**NOTE:** Read this manual carefully. It contains important information about safety during installation, use and maintenance. If the equipment is utilized in a manner not indicated by the manufacturer, the safety provided by the equipment may be adversely affected.

## Meaning of the symbols used in relation to the machine

	Danger. Danger of injuries.
	Read these instructions in their entirety before putting the device into operation.
	Risk of electric shock.
	Caution: the surface of the lid and the pot, the water and the steam at 100 °C, can cause burns. Use suitable protection.
	Water outlet connection.
	Pressurised air inlet connection.
	Outlet of steam at 100 °C / 212 °F.

## Precautions

- Keep this manual in a suitable place for future reference.
- Check the general conditions of the device immediately after removing the packaging.
- Breakdown caused by incorrect installation or incorrect use of the device is excluded from the guarantee.
- To avoid the risk of electric shock, this device must only be connected to a power supply with proper earthing. Before connecting the device, check that the identification label indicates the correct voltage. The installation must comply with all applicable safety standards.
- Remove any remains of dirt or water from the connection of the cable.


- Avoid excessive traction, twisting or folding of the cable for connection to the electricity grid.
- Do not touch the device or the cable with wet or moist hands.
- Install the product away from sources of heat, dust, high humidity, and from splashing of water or inflammable material (alcohol, solvents, etc.).
- Unplug the device immediately if you detect noise, smells or smoke.
- Ensure that children and unspecialized staff do not handle the device.
- If in doubt, do not put the device into operation. Contact a qualified and authorized repairer.
- Do not disassemble, repair or modify the device.
- Unplug the device if there is a storm, and when out of use for a long time.
- Unplug the device before cleaning it, do not wet it directly, and do not use chemical products.
- Contact your distributor if you have any questions or for detailed information about problem solving.
- Before using the device, ensure that it is positioned on a firm and solid base.
- The cable can only be replaced by one with identical characteristics to that supplied by the manufacturer.
- Protect yourself properly; use gloves and safety glasses.
- Do not allow anybody to bridge the safety systems.
- Do not use the machine for any use other than its intended purpose.

## Intended use


The Ref. 030425 is a polymerisation pot, easy to handle and especially designed for the polymerisation of prostheses and other elements within the dental prosthesis sector.

## Installation

1. Carefully unpack the unit. Check that it has not suffered any deterioration during transport.
2. Place the device on a firm, flat surface away from any focus of heat or vibrations.
3. Connect one end of the Ø 4x6 tube to the air inlet (H). Connect the other end to a source of pressurised air (minimum 3 bar).
4. Choose a position for the drainage bottle in a place located at least 20 cm below the bottom of the machine. The water expelled from the tank is at boiling temperature, which means that an accidental fall of the bottle can cause an accident. Therefore, take care that the position of the drainage bottle ensures optimum stability and safety. Once the bottle is in position, connect the Ø 6x8 drainage tube to the outlet (G), and insert the other end into the upper orifice of the bottle.

- An accessory is available to work at 1.5 bar (L).  
Ref. 030425-45.
- To discharge water, open the water valve (F) situated on the right side of the machine (as long as the bottle has been correctly connected).

- Connect the device to a suitable electrical connection with proper earthing.

## Use of the keys

Program selection.

Edition of time and temperature values.

Access to edition of temperatures.

When the cycle is underway, holding the button down for 3 seconds enables you to visualise the current temperature.

Access to edition of times.

When the cycle is underway, holding the button down for 3 seconds enables you to visualise the time remaining for the phase underway.

Short push: cycle startup.

Long push: end of cycle.

Activate and deactivate the entry of pressure.

## Polymerisation programs


There are three types of program: short, long and continuous.

### Short programs:

Consist of a consigned polymerisation temperature  $^{\circ}\text{C1}$  and a maintenance time  $t^{\circ}\text{C1}$ .

### Long programs:

Consist of a preheating temperature  $^{\circ}\text{C1}$ , a maintenance time  $t^{\circ}\text{C1}$ , a temperature rise slope  $r^{\circ}\text{C2}$ , a consigned polymerisation temperature  $^{\circ}\text{C2}$  and a maintenance time  $t^{\circ}\text{C2}$ .

Short program: if  $^{\circ}\text{C2}$  is programmed lower than  $^{\circ}\text{C1}$ .

Long program: if  $^{\circ}\text{C2}$  is programmed higher than  $^{\circ}\text{C1}$ .

### Continuous program:

This is a short program which keeps temperature in order to let any user starting a new polymerisation cycle at any time.

## Default programs

The pot leaves the factory with the following programs:

N°	Pressure	$^{\circ}\text{C1}$	$t^{\circ}\text{C1}$	$r^{\circ}\text{C2}$	$^{\circ}\text{C2}$	$t^{\circ}\text{C2}$
P1	Auto	45	0h15	-	-	-
P2	Auto	55	0h15	-	-	-
P3	Manual	100	0h30	-	-	-
P4-L	Manual	70	1h30	0h30	100	0h30
Pc	Auto	55	0h15	-	-	-

$Pc$  This is the continuous program (see in its section).

$P3$  Means that the program P3 is short (one temperature sign).

$P4-L$  Means that the program P4 is long (two temperature signs).

## Program selection

- Press the main switch (K) and wait until the display stabilises.

P1 Program selected (short).

45° Consigned temperature.

P4-L Program selected (long).

70°C Consigned temperature  $^{\circ}\text{C1}$ .

100°C Consigned temperature  $^{\circ}\text{C2}$ .

4. Use keys  and  to select one of the 5 programs (Pc, P1, P2, P3 or P4-L).
- Pc* This is the continuous program (see in its section).
- P3* Means that the program P3 is short (one temperature sign).
- P4-L* Means that the program P4 is long (two temperature signs).

## Modification of temperatures

**Note:** Programs cannot be edited whilst underway.

1. Select the program you wish to modify.
2. Press the key  to modify the temperature °C1 of the program selected.

°C 1 70°	Consigned temperature.
-------------	------------------------

3. Press the key  again to modify the second temperature °C2.

°C 2 100°	Consigned second temperature.
--------------	-------------------------------

4. Use the keys  and  to modify the values.
5. Press the key  again to return to program selection. The values modified will be saved in the memory automatically.

**Note:** If you have selected a short program and you wish to program a long one, hold down the key  whilst °C1 is showing in order to gain access to °C2. Program a °C2 temperature higher than °C1, and you will gain access to the rest of the parameters for the long program.

## Modification of times

**Note:** Programs cannot be edited whilst underway.

1. Select the program you wish to modify.
2. Press the key  to modify the maintenance time t °C1.

t °C 1 1h30	Maintenance time for the first temperature.
----------------	---

3. Use the keys  and  to modify time t °C1.

4. Press the key  again to modify the rise slope to the second temperature r °C2.

r °C 2 0h30	Time for temperature rise slope to the second temperature.
----------------	--

5. Press key  again to modify the maintenance time t °C2.

t °C 2 0h30	Maintenance time for the second temperature.
----------------	--

6. Use the keys  and  to modify the values.
7. Press the key  again to return to program selection. The values modified will be saved in the memory automatically.

**Note:** If you have selected a short program and you wish to program a long one, hold down the key  whilst °C1 is showing in order to gain access to °C2. Program a °C2 temperature higher than °C1, and you will gain access to the rest of the parameters for the long program.

## Operation

1. Open the lid and pour water into the tank.

	To open the lid, press the red button on the handle and turn to the right.
--	--

2. Activate the main switch located on the back of the device (K). The green light will come on.
3. Select the program you wish to run.
4. Press key . The program will start up.

The pilot light next to key  will blink, indicating that it is heating up, and will remain with a constant light once the consigned temperature is reached.

5. To terminate a program which is underway, hold down the key  for a few seconds.

	When the cycle is underway, holding the key down for 3 seconds enables you to visualise the current temperature.
	When the cycle is running, holding the key down for 3 seconds enables you to visualise the time remaining for the phase underway.

6. When the assigned temperature is reached, a beep will sound and the machine will show the message **L<sub>o</sub>Ad.**

Insert the material for polymerisation and close the lid of the pot.


Caution: the surface of the lid and tank, the water and the steam at 100 °C can cause burns. Use suitable protection.

Press  to continue the cycle.

7. When the polymerisation temperature is reached, proceed to manual or automatic pressurisation (see the following section “*Pressurisation*”).
8. When the time is complete, the depressurisation phase will begin and the message **dEP<sub>r</sub>** will be shown.

**Note:** To discharge water, open the water valve (F) situated on the right side of the machine (as long as the bottle has been correctly connected).

9. After 2 minutes of depressurisation, a beep will sound and the machine will show the message **End.**
10. Before opening the lid, check that the water discharge is fully complete and that the pressure indicated on the pressure gauge is 0.

Caution: the surface of the lid and tank, the water and the steam at 100 °C can cause burns. Use suitable protection.

11. Press  to erase the message **End.**
12. Once the job is complete, turn off the switch (K).

## Pressurisation

The unit is equipped to regulate the pressure to approximately 2 bars (it does not matter if this varies slightly).

**Note:** A small outlet of air from the covering of the lid can be observed. This is completely normal. This is from a valve which regulates the pressure inside the tank.

### Automatic pressurisation

If the consigned polymerisation temperature is equal to or less than 60 °C, pressurisation will be automatic (as long as the temperature has been reached and the user has indicated that the item for work has been inserted).

### Manual pressurisation

If the polymerisation temperature is higher than 60 °C, the user will have to carry out the pressurisation manually.

At the suitable moment, a beep will sound and the machine will show the message **bAr.**

The user will have to press the key  to introduce the pressurised air into the tank.

Do not pressurise the tank with the lid open.

## Continuous program operation (Pc)

This is a short program which keeps temperature in order to let any user starting a new polymerisation cycle at any time.

**Note:** The maximum temperature allowed in the Pc program is limited to 75 °C to reduce the possibility of the tank running out of water due to evaporation.

### Limit time of continuous operation

Since the Pc program starts, the temperature will be maintained for a time limit of 10 hours. If you wish, you can modify this value between 0 and 15 hours:

- Select the program Pc.
- Press the key  to modify the maintenance time t°C1.

t°C1  
0h 15

Maintenance time for the first temperature.

- Now, keep pressing the key  to modify the limit time of continuous operation tPc.

tPc  
10

Limit time of continuous operation (hours).

- Use the keys  and  to modify values.
- Press the key  again to return to program selection. The values modified will be saved in the memory automatically.

## Cycle end warning – °C/°F

By default, the machine will emit an audible warning for **5 seconds** at the end of the polymerization.

- Simultaneously press keys  and  for 2 seconds.

**6EEP**  
**5** Cycle end warning time (seconds).

2. Use the keys  and  to modify values (0-30 sec). If you select the value "ON", warning will remain indefinitely until a key is pressed.

3. Press again the key  to validate the value and move on to the next parameter.

**dE9r**  
**°C** Temperature in Celsius degrees.

4. Use the keys  and  to change between °C and °F.

5. Press again the key  to validate the value.

## Error messages

<b>E-01</b>	Broken interior thermal probe.	*
<b>E-02</b>	Broken resistor sensor.	*
<b>E-10</b>	The safety thermostatic switch actuated.	*
<b>E-20</b>	Interior overheating	*
<b>E-21</b>	Resistor overheating	*
<b>E-22</b>	Electronic controls overheating	*
<b>E-80</b>	Electricity supply cut	**
<b>E-81</b>	Lack of water	**
<b>E-82</b>	Excessive waiting time	**

\* Technical service actuation is needed.

\*\* Press any key to erase the error message.

## Safety elements

The machine has the following safety elements:

- Air inlet pressure regulator.
- Safety valve on lid.
- Interior safety valve.
- Detection of anomalies by electronic control.
- Resistor temperature sensor.
- Safety thermostat.

## Maintenance and cleaning

To ensure optimum reliability and long life of the device, we recommend:

- After each use, we recommend cleaning away remains of materials from inside the tank with normal washing up liquid. In order to do so, remove the plastic plate from the bottom of the tank and clean the inside of the tank, taking care not to damage the temperature probe. Rinse the inside of the tank with abundant water.
- To prevent remains of material from obstructing the exit from the tank, clean the inside of the tank periodically (approx. every 5 uses). To do so, heat the water to a temperature of 100 °C. Then empty the tank through to the drainage bottle.
- Periodically, and depending on water hardness, clean away remains of limescale which may have been deposited on the bottom of the tank.
- Do not use solvents or any other inflammable product.
- Proceed with caution when handling the lid in order not to bash it (there are safety valves inside).
- Replace the rubber seal of the lid periodically. For this operation, we recommend using original MESTRA replacements.
- If you feel excessive resistance when turning the lids, or it does not seal well, you can clean the rubber seal and lubricate it lightly with talcum powder.

## Technical characteristics

Height: ..... 35 cm  
 Width: ..... 36 cm  
 Length: ..... 40 cm  
 Weight: ..... 5,5 kg  
 Pressure: ..... 2 bar  
 Maximum temperature: ..... 100 °C  
 Tank capacity: ..... 6 L  
 Tank dimensions: ..... Ø22 x 15 cm  
 Heating power: ..... 600 W

## UE statement of conformity


The company **Talleres Mestraitua S.L.**  
 with social address: Txori-erri Etorbidea, 60  
 48150 Sondika - Bilbao - Spain

hereby declares that the product

**Automatic polymerisation pot 030425**

Meets the requirements set out in the following directives:

- 2014/35/UE (low voltage)
- 2014/30/UE (electro magnetic compatibility)
- 2006/42/CE (machinery safety)

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