CEN-MOTOR SM-2

R-040064

ESPAÑOL	4
ENGLISH	7
FRANÇAIS	10
DEUTSCH	13
PORTUGUÊS	16
ITALIANO	19
NORSK	22

MESTRA®

TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

T +(34)944530388 - F +(34)944711725


E-mail: mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA


Rev. 09/01/13

A	Brazo / Arm / Bras / Arm / Braço / Braccio / Arm
B	Tapa superior / Lid / Couvercle supérieur / Oberer Deckel / Tampa superior / Coperchio superiore / Toppdeksel
C	Interruptor / On switch / Interrupteur / Schalter / Interruptor / Interruttore / Hovedbryter
D	Control / Control / Contrôle / Steuerung / Controlo / Controllo / Styreenhet

- | | |
|----|--|
| 1. | Cilindro / Cylinder / Cylindre / Zylinder / Cilindro / Cilindro / Støpesylinder |
| 2. | Crisol / Crucible / Creuset / Tiegel / Cadinho / Crogiolo/ Digel |
| 3. | Tuerca de fijación / Fixing nut / Ecrou de fixation / Befestigungsmutter / Porca de fixação / Dado di fissaggio / Festemutter |
| 4. | Guía / Guide / Guide / Gegengewicht / Corrediza / Guida / Motvekt |
| 5. | Contrapeso / Counterweight / Contrepoids / Führung / Contrapeso / Contrappeso / Føringsskinne |

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The CEN-MOTOR SM-2 is a modern centrifuge especially designed for centrifuge moulding small metal parts for the dental prothesis sector. To ensure optimum operating conditions and long life for the unit, please read this manual carefully.

INSTALLING THE UNIT

1. Carefully unpack the CEN-MOTOR SM.-2 and check that it has not suffered any damage in transit.
2. Set it up on a firm, flat surface away from any source of vibration or heat.
3. Set the unit up on its four supporting feet, then adjust the position of those feet via the adjusting nuts. To work properly the CEN-MOTOR SM-2 must be perfectly horizontal, so a spirit-level should be used in setting it up.
4. Connect the CEN-MOTOR SM-2 to an earthed 230 V, 50/60 Hz power outlet.

OPERATION

- ? Connect the unit to a 230 V, 50/60 Hz earthed socket.
- ? Fit the cylinder (1) into its housing, ensuring that the outlet is not obstructed and is perfectly aligned with the hole in the crucible. Do this by using the lifting device, which has a scale matching the diameter of the crucible. This simplifies the job of fitting the cylinder onto the crucible. We recommend the use of MESTRA cylinders and crucibles, which are designed to give maximum performance with the CEN-MOTOR SM-2.
- ? Fit the arm counterweight (4) at a suitable position. This counterweight slides easily along the guide (5) when the fixing clip is released. To facilitate this task, the guide is marked with numbers that match the diameter of the cylinder. (Read the table guide at the end of this manual).
- ? The arm has an automatic device to retract the counterweight and cylinder/crucible assembly around two articulations placed in the center of the unit. Please, check out that this device have worked right. Incorrect retraction (i.e. failure to articulate fully) could lead to the loss of the pour. Also ensure that the articulations can turn freely. If any resistance is encored, consult your authorized **MESTRA** dealer immediately.
- ? Retract the counterweight and cylinder/crucible assembly by turning the two articulations

towards the center of the unit, so that their final shape is as shown in the figure below:

- ? Once the arm is counterweighted, remove the cylinder and crucible from their housing and pre-heat them in the oven (to prevent heat shock). We recommend MESTRA series HP ovens.
- ? Wear suitable protective goggles and gauntlets when handling blowtorches, cylinders, hot crucibles and the oven. After pre-heating fit the crucible into its housing on the arm.
- ? Place the metal pellets to be melted in the crucible. Do not use too much: this could lead to crucible overflow and molten metal splashing. To prevent heat-shock problems, pre-heat the crucible before melting pellets.
- ? Melt the metal in the crucible with an oxyacetylene torch. For optimum results continue to apply heat to the molten metal for a further 10 seconds for normal melting point precious metals and 15 for those with high melting points after melting is completed.
- ? When the melt reaches the right temperature, take the pre-heated cylinder from the oven and place it in its housing on the arm, ensuring that it is properly aligned and attached.
- ? With the cover open, use the on switch (6) on the CEN-MOTOR SM-2.
- ? The display (12) shows the centrifuge time. Use the up (10) and down (11) keys to select the time (the unit allows selection between 10 and 60 seconds only).
- ? Press the TORQUE (8) key to select one of the possible delayed start-ups. When the key is pressed the delayed start-up selection changes and the LED lights up. If "0" start-up is selected, the motor will quickly increase to 100 % of its speed; if "4" start-up is selected, the motor will slowly increase to 100 % of its speed.
- ? At the end of the manual there is a table with recommended TORQUE values, in accordance with the diameter of the cylinder and the level of smelting.
- ? Finally, to memorize the option selected (start-up type and time) and commence the centrifuge cycle, press the START key (7). The pilot will light up and the centrifuge is ready for operation.

- ? Close the cover. The CEN-MOTOR SM-2 will start up and fill the cylinder with molten metal.
- ? Wait around 1 minute for the metal to solidify in the cylinder.
- ? Lift the lid. The unit will stop immediately, and the cylinder and crucible can be removed.

MAINTENANCE AND CLEANING

To ensure reliable operation and a long lifetime for your unit, proceed as follows:

- ? Brush out the inside of the housing after each use to remove any splashes of metal.
- ? Remove any metal splashes on the crucible-holder arm, ensuring particularly that the articulations turn freely. Grease the articulations from time to time with mineral grease.
- ? To clean the housing of the unit use a cloth moistened in soapy water.
- ? Never use alcohol or other inflammable liquids to clean the CEN-MOTOR SM-2.
- ? It is normal for the crucible-holder arm and drive shaft to oscillate while the CEN-MOTOR SM-2 is working. Rocking movements of up to 15 mm are acceptable and envisaged by **TALLERES MESTRAITUA, S.L.** If bigger movements occur or unusual noises are heard, consult your **MESTRA** dealer.

PRECAUTIONS

- ◆ Before connecting the unit ensure that the power supply is 230 V, 50/60 Hz and earthed.
- ◆ Do not allow children or unskilled personnel to handle the unit.
- ◆ Do not over-fill the crucible: excess molten metal may splash out.
- ◆ Use suitable goggles and safety gauntlets when melting metal and handling cylinders and crucibles.
- ◆ Ensure that the fixing nut on the drive shaft is tight before starting up the unit.
- ◆ Before using the unit ensure that it is properly set up on a solid, level base.

- ◆ Check that the crucible outlet and cylinder runner are properly lined up and attached.
- ◆ Ensure that the crucible holder arm is properly counterweighted and folded before starting up the unit.

- ◆ The machine should spin anticlockwise.
- ◆ Do not place liquids or inflammable materials in the crucible.
- ◆ If the motor shaft locks during normal operation, disconnect the unit from the mains before attempting to release it by turning it manually.
- ◆ After the pour, handle the cylinder and crucible with care: they are hot enough to burn.

HOW TO FIT THE COUNTERWEIGHT

Here is a table with the position of the counterweight depending on the cylinder MESTRA that you are using. When graduating the counterweight, make sure that the line located just above the correspondent inscription on the guide is covered. For example, if a 70 mm cylinder is placed, you should see the word marking 70-9X and the line above this inscription should remain covered.

CILINDER	INSCRIPTION
30 mm	30-1X
40 mm	40-3X
50 mm	50
60 mm	60-6X
70 mm	70-9X
80 mm	80-CEG
1 X	30-1X
3 X	40-3X
6 X	60-6X
9 X	70-9X
Metal Protheses little	CEP
Metal Protheses medium	CEM
Metal Protheses big	80-CEM

TORQUE GRADUATION TABLE			
Cylinder diameter	Low melting point metals	Medium melting point	High melting point
30	4	4	3
40	4	4	3
50	3	3	2
60	3	2	1
70	1	1	0
80	1	1	0

TECHNICAL CHARACTERISTICS

Single-phase motor:	370 W
Service consumption:	230 V, 50/60 Hz
Height:	320 mm
Width:	600 mm
Depth:	560 mm
Weight:	40 kg

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