## R-080010 / R-080020 / R-080025

E.....	ELECTROVIBRADORES.....	4
GB.....	ELECTRO-VIBRATORS.....	5
F.....	ELECTRO-VIBRATEURS.....	6
P.....	ELECTRO-VIBRADORES.....	7
I.....	ELETTROVIBRATORE.....	8
D.....	ELEKTRISCHE RÜTTLER (VIBRATOREN).....	9

**MESTRA®**



TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60


Tel. +34 944530388 - Fax +34 944711725

E-mail: [mestra@mestra.es](mailto:mestra@mestra.es) - [www.mestra.es](http://www.mestra.es)

48150 SONDIKA - BILBAO - ESPAÑA

R-080010

R-080020 / R-080025

- A Ventosas de base / **Base suction caps** / *Ventouses* / Ventosas de base / **Ventose base** / *Saugfüße*
- B Regulador de vibraciones / **Vibration regulator** / *Régulateur de vibrations* / Regulador de vibrações / **Regolatore vibrazioni** / *Schwingungsregler*
- C Plataforma de vibración / **Vibration platform** / *Plateau de vibration* / Plataforma de vibração / **Piattaforma di vibrazione** / *Rütteltisch*
- D Interruptor general / **On-off Switch** / *Interrupteur Marche-Arrêt* / Interruptor geral / **Interruttore generale** / *Hauptschalter*

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The R-080010 / R-080020 / R-080025 is a practical electro-vibrator unit specially designed for compacting materials for the jewellery and dental mechanics sectors. To obtain best performance from the unit, we recommend you read the following instructions carefully.

- Place the unit on a smooth, firm, perfectly horizontal surface. The electro-vibrator is fitted with rubber suction caps (A) to ensure firm attachment.
- Connect the plug to a 220 V, 50/60 Hz earthed power outlet.
- Place the cylinder or items to be vibrated on the unit's vibration platform. In the case of impression bowls, the mixture can also be compacted by resting the bowl firmly on the edge of the unit and holding down the bowl handle with your hand.
- Adjust the vibration regulator (B) to its minimum setting (turn the knob anti-clockwise), then connect the On switch (D) and the electro-vibrator will start up.
- Slowly move the unit's regulator (B) clockwise and check that the oscillation rate gradually increases. Set the regulator to a suitable amplitude for the operation.
- After a short time, once the mixture is properly compacted, switch the unit off using the On/Off switch (D).

## MAINTENANCE AND CLEANING

To ensure a long, trouble-free service life and performance of the unit, we recommend as follows:

- After using the unit, always clean away any splashes of mixture before they dry.
- Clean the vibration platform and outer housing of the electro-vibrator regularly with a damp cloth. The rubber cover protecting the platform can be easily removed to facilitate cleaning.
- As a precaution, always disconnect the unit from the mains before cleaning it.
- To remove the vibration platform (C), remove the 4 bolts next to the suction pads (A) on the base of the unit (see photograph, page 2).

## PRECAUTIONS

- ◆ Before connecting the unit ensure that the power supply is 220 V, 50/60 Hz and earthed.
- ◆ Do not use the unit in wet places or close to heat sources.
- ◆ Do not place easily broken containers or other objects on the vibration platform.
- ◆ Do not allow children or unskilled personnel to handle the unit.
- ◆ Do not pull, twist or bend the mains cable excessively.
- ◆ Before cleaning the unit, ensure that it is unplugged.
- ◆ Never immerse the unit in liquid when cleaning it.
- ◆ Do not place hot containers or objects on the vibration platform, as this may damage the rubber.

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