R-080251

E.....	MÓDULO DE MICROCHORREADO.....	3
GB.....	MICRO-BLASTING MODULE.....	6
F.....	MODULE DE MICRO-SABLAGE.....	8
P.....	MÓDULO DE MICRO JACTO.....	10
I.....	MODULI DI MICROSABBIATURA.....	13

MESTRA[®]


TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tfno. + 34 94 453 03 88 Fax + 34 94 471 17 25

E-mail: mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

- | | |
|---|--|
| A | Tubo de conexión a red de aire / Compressed air connection tube / <i>Tube raccord à la prise d'air</i> / Tubo de conexão a rede de ar / Tubo di collegamento rete d'aria |
| B | Manorregulador / Pressure regulator / <i>Mano-régulateur</i> / Mano regulador / Regolatore |
| C | Manómetro / Pressure gauge / <i>Manomètre</i> / Manómetro / Manometro |
| D | Pedal de control / Control pedal / <i>Pédale de fonction</i> / Pedal de controlo / Pedale di controllo |
| E | Tapa / Tank stopper / <i>Bouchon</i> / Tampa / Coperchio |
| F | Depósito de esmeril / Grinding medium tank / <i>Silo</i> / Depósito de esmeril / Serbatoio di smeriglio |

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The R-08251 model is a practical sandblasting module, specially designed to microblast small components of dental prostheses.

The unit should be used with a microblasting booth to prevent the dispersion of dust and grinding particles thrown out during operation.

INSTALLATION

- Carefully unpack the unit.
- Micro-blasting modules should be used with booths equipped with exhausters. Set up the unit on a flat, rigid surface away from any sources of vibration or heat, and in an environment which is not too humid. Too much humidity could lead grinding particles to compact in the booth and clog the micro-blasting nozzles.
- Connect the tube provided (A) to an air outlet at a pressure of no more than 5 bar.
- Locate the pedal of the micro-blasting gun (D) within reach of your foot, in a position where it can be worked comfortably.
- Remove the threaded stopper on the tank (E) and fill the vessel with a grinding medium with a particle diameter of less than 250 microns. The grinding medium should be at least 30 mm below the rim of the brass tube. Grinding medium heights of less than 35 mm are not recommended (see figure). Screw the tank stopper back in place, ensuring that the seal is leak-tight.

PRESSURE REGULATION

The blowing pressure of the micro-blasting gun is the same as the pressure at the air outlet (compressor). This pressure may be regulated between 2 and 5 bars to suit the nature and characteristics of the work in hand.

PRESSURE SHOULD NEVER EXCEED 5 BARS. To adjust the pressure use the pressure regulator, which limits the output pressure from the compressed air unit.

To adjust the pressure in the micro-blasting module (R-080250 only), pull up the pressure regulator knob (B) beside the grinding medium vessel and turn it to the desired setting.

The pressure is shown on the pressure gauge (C) built into the regulator. Once pressure is selected, press the regulating knob back down to lock it.

USING THE UNIT

To micro-blast workpieces correctly, proceed as follows:

1. Check that the level of grinding medium in the micro-blasting module vessel is within the limits for work.
2. Use suitable protective goggles and a face-mask, and work in a booth equipped with an exhauster.
3. Start up the exhauster.
4. To use the micro-blasting gun, press the operating pedal.

5. Point the sand jet towards the workpiece. Do not allow it to blast objects other than those to be treated.
6. When work is completed, release the pedal. Keep the exhauster working for a few seconds to remove all dust particles from inside the booth, then stop it.

MAINTENANCE AND CLEANING

To ensure a long service life for the unit, proceed as follows:

- Always use high-quality grinding medium. **Talleres Mestraitua, S.L.** recommends **MESTRA** medium. The gun on the micro-blasting module is factory supplied with a nozzle diameter of 1.25 mm and can handle grinding medium up to 250 microns in diameter.
- Use dry piston air compressors whenever possible. The MESTRA R-080340 and R-080350 compressors are suitable for this task. Oil particles in the blowing air can cause lumps of sand which will clog the gun nozzles.
- Do not fill the micro-blasting module vessel with the grinding medium from the booth floor: large diameter particles (impurities appearing during workpiece blasting) may clog the nozzle of the gun.
- For general cleaning of the unit, use a cloth moistened in soapy water.

If the gun becomes clogged, dismount the nozzle-holder and remove the nozzle from inside. Unclog the nozzle with a pin or fine wire and re-mount it on the gun.

The micro-blasting gun is factory supplied with a 1.25 mm diameter nozzle. However, 1 mm and 0.75 mm nozzles can also be fitted: consult your **MESTRA** dealer. The following table gives grinding medium particle diameters for the different nozzle diameters.

MAXIMUM PARTICLE DIAMETER FOR EACH NOZZLE	
Nozzle	Particle diameter
1.25 mm	250 microns
1 mm	150 microns
0.75 mm	75 microns

PRECAUTIONS

- ◆ Never exceed a pressure of 5 bars.
- ◆ Always work with a micro-blasting booth, and use protective gauntlets, goggles and a face-mask.
- ◆ Use an exhauster while micro-blasting.
- ◆ Do not allow children or unskilled personnel to handle the unit.
- ◆ Use only guaranteed high-quality grinding media (we recommend MESTRA media). Grinding media can be toxic if dust is inhaled.
- ◆ Ensure that the compressed air system used to power the unit is properly protected by at least one pneumatic discharge valve.
- ◆ If leaks are detected in the compressed air system, consult your installer immediately.
- ◆ After filling the module's tank with grinding medium, ensure that it is tightly sealed.

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