
R-080262 / R-080263 / R-080264
R-200262 / R-200263 / R-200264

E.....	CHORREADORAS DE ARENA CONSTELLATION	3
GB ...	CONSTELLATION SANDBLASTERS	6
F.....	SABLEUSES CONSTELLATION	8
P.....	JACTO DE AREIA CONSTELLATION	11
I.....	SABBIATRICI CONSTELLATION.....	14
D.....	CONSTELLATION SANDSTRAHLER	19

MESTRA®

TALLERES MESTRAITUA S.L.
 Txori-Erri Etorbidea, 60
 48150 SONDIKA - BILBAO - ESPAÑA
 Tel. +34 944530388 - Fax +34 944711725
 mestra@mestra.es - www.mestra.es

Rev. 27/12/12

A	Vaso depósito de esmeril / Emery pan / <i>Vase réservoir d'alumine</i> / Vaso depósito de esmeril / Vaso serbatoio smeriglio / <i>Schmirgelbehältergefäß</i>
B	Manorregulador vaso / Pan pressure adjustment knob / <i>Mano régulateur de réservoir</i> / Regulador de pressão vaso / Regolatore pressione vaso / <i>Gefäßdruckregler</i>
C	Manómetro de presión / Pressure gauge / <i>Manomètre de pression</i> / Manómetro de pressão / Manometro pressione / <i>Druckmanometer</i>
D	Interruptor luz cabina / Chamber light switch / <i>Interrupteur de la lumière de la cabine</i> / Interruptor luz cabina / Interruttore luce cabina / <i>Kabinenlichtschalter</i>
E	Conmutador de selección de vaso / Pan selection switch / <i>Commutateur de sélection de vase</i> / Comutador de selecção do vaso / Commutatore di selezione vaso / <i>Gefäßwahlumschalter</i>
F	Manguitos de tela / Sheaths / <i>Manchons en tissus</i> / Mangas de tela / Manicotti di tela / <i>Gewebeschläuche</i>
G	Tubo de entrada aspiración exterior / External suction inlet tube / <i>Tube pour aspiration extérieure</i> / Tubo de entrada aspiração exterior / Tubo d'entrata aspirazione esterna / <i>Einlassrohr äußere Absaugung</i>
H	Aire comprimido / Compressed air / <i>Air comprimé</i> / Ar comprimido / Aria compressa / <i>Druckluft</i>

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The Constellation sandblaster family is composed of several sandblasting chambers specially designed for precision-blasting small parts used in the dental prosthesis industry. For optimum output and performance of this unit, we advise you to read the following instructions with great care:

INSTALLING THE UNIT

- Carefully unpack the unit. You will find the following fittings inside the chamber:
 - An air connection tube, diameter 6 mm.
 - A 90 ° PVC elbow fitting to connect the outgoing tube to suction if required.
 - Pedal drive mechanism.
- A rigid flat surface must be chosen for installation of the unit, at a safe distance from any sources of vibrations or heat, and the surrounding atmosphere must not contain excess moisture, since this could cause compacting of the emery particles inside the chamber, and lead to obstructions in the microblasting nozzles. The unit must also be positioned on a level which allows objects inside the chamber to be handled with ease.
- Hook up the air connection tube supplied with the unit to the connection at the rear of the unit. The other end must be fitted to a compressed air facility, with at least 7 bars of pressure and minimum flow rate 100 litres/minute.
- Hook up the electric plug to 230 V, 50/60 Hz, which must be earthed.
- Connect the pedal drive tube to the coupling nozzle on the right-hand side of the unit. Then position the pedal for comfortable foot operation. Depending on operational requirements and/or handling comfort, it is possible to position the pedal for knee operation.
- Unscrew the sandblasting module covers (A), and fill the pans with emery. During this operation, **it is most important not to fill the pans to more than ¾ of their capacity.** Put the covers back on, and ensure they have been screwed on tightly.
- In relation to the calibre of the nozzles, the machine is supplied ex-works to the following criterion of nozzle diameters and colours:

Yellow nozzle	Ø 0.75 mm	
Blue nozzle	Ø 1 mm	
Green nozzle	Ø 1 mm +	
Red nozzle	Ø 1.25 mm	

REGULATING THE PRESSURE

The blowing pressure of the blasting guns can be regulated depending on the nature and characteristics of the operation to be carried out. Use the pressure-adjustment knob (B) on the top right of the unit – pull the knob upwards, and rotate it to the required pressure. The unit has a pressure gauge (C) for easy selection of operating pressure. When the pressure has been adjusted, lock out the regulator knob by pushing it down to the initial position.

Pressure gauge (C) shows the blasting module gun blow pressure at all times, even when the gun is non-operational.

CONNECTION TO THE SUCTION UNIT

The microblasting chamber must be permanently connected to forced suction via an outgoing tube on the top left of the unit (G). Suction must be able to produce a slight depression inside the chamber for continuous gentle circulation of the powder generated towards the outgoing tube.

Depending on the requirements and geometry of the external suction unit, it is possible to use the PVC elbow fitting supplied for specific needs.

USING THE UNIT

The following operations are recommended for proper precision-blasting:

- Open the lid of the chamber, and put in the part you wish to blast. Parts may also be placed inside the chamber via sleeve (F). Check there is enough emery in the pans.
- Close the lid of the chamber, and ensure it is tightly shut (the glass should fit snug against the housing).
- Suitable protection gear must be used: goggles, gauntlets and a mask.
- Start up the forced suction facility. The light in the chamber may be activated for a better view of the part, if required, using switch (D).
- Use the main switch (E) to select the pan. The colour coding system used makes it easy to see which is the active pan at all times.
- Pick up the mobile gun selected inside the chamber. The colour of the gun's nozzle-holder matches the colour on the pan covers (red, blue, yellow, green). When the gun has been selected, press the drive pedal to begin sandblasting.
- Direct the blast carefully to prevent the stream hitting the sides of the chamber directly, since this could damage the

chamber. If you wish to test the blasting pressure or strength, direct the stream at the rubber dolly inside the chamber.

8. After the blasting operation, release the pedal, stopping the sand flow and activating immediate air discharge to relieve pressure in the pans. Then switch off the suction unit, open the lid of the chamber, and take out the part.

MAINTENANCE AND CLEANING

To ensure a lengthy service life, we recommend the following:

- Always use approved quality emery - **Talleres Mestraitua, S.L.** recommends **MESTRA** emery. The mobile blasting guns with nozzle diameters 1.25 mm and 1 mm can respectively blow emery of maximum diameter 250 and 150 microns. The mobile blasting gun with nozzle diameter 0.75 mm can blow a maximum diameter of 75 microns (see the table).
- It is preferable to use dry-piston air compressor units, since oil in the air could cause the sand to compact and obstruct the gun nozzles.
- Never fill the modules with emery deposited on the floor of the chamber, since particles with large diameters could obstruct the nozzles, and this emery is now contaminated with particles eroded off the parts blasted.
- The silencer filters inside the pans must be replaced every 5,000 or 10,000 operations (around 1 or 2 years). The incoming filter fittings on each pan must be replaced every 500 or 1,000 hours (much depends on the quality and metrage of the sand used). Original MESTRA spare parts must always be used.
- To prevent any abrasion wear on the walls inside the chamber, take care not to direct the stream of emery directly onto the walls.
- Clean out the emery inside the chamber on a periodic basis, using the practical emptying plug inside.
- Any water which has condensed inside the incoming air filter must be bled out regularly. The chamber has a window to view the amount of water in the filter.
- Use a dry cloth or a cloth moistened with soapy water for general cleaning of the unit.
- If the guns become obstructed, remove the nozzle-holder and take out the nozzle. Unblock the nozzle by blowing in the opposite direction to the flow of air, or use a pin or a thin wire to unblock it, and then refit the attachment.

- The blasting module guns on the 2-pan Constellation model are supplied ex-works with nozzles 1 and 1.25 mm in diameter. 0.75 nozzles, however, may also be used with this model (consult your **MESTRA** distributor).
- The maximum diameter of emery particles in relation to nozzle diameter is shown in the table below:

MAXIMUM PARTICLE DIAMETER / NOZZLE	
Nozzle	Diameter
1.25 mm	Up to 250 microns
1 mm	Up to 150 microns
0.75 mm	Up to 75 microns

PRECAUTIONARY MEASURES

- ♦ Never exceed 5 bars of pressure. Optimum working pressure is between 2 and 4 bars.
- ♦ Operations should always be carried out with the lid of the chamber closed, with suitable protection. Masks and goggles should be worn.
- ♦ Connect the chamber's outgoing tube to a suction facility.
- ♦ The protection sleeves (F) must be replaced immediately following any tearing or damage.
- ♦ Do not allow children or unqualified staff to operate the unit.
- ♦ Use only guaranteed-quality materials as emery (we recommend **MESTRA** emery). Poor quality emery may be toxic when inhaled, and may also cause obstructions in the blow nozzles.
- ♦ Hook up the power supply cable to 230 V, 50/60 Hz, earthed.
- ♦ Consult your fitter immediately if any leakages are observed in the compressed air circuit.
- ♦ Before carrying out any operation, first ensure that the micro blasting module covers are closed tightly.

TECHNICAL CHARACTERISTICS

Height:	320 mm
Length:	390 - 510 mm
Width:	440 mm
Weight:	12.5 - 14 - 15.5 kg
Working pressure:	2-6 bar
Voltage:	230 V, 50/60 Hz
Power consumption:	4 W

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