R-080400 / R-080402

E.....	POLIMERIZADORA	4
GB.....	POLYMERISER	9
F.....	POLYMÉRISEUSE	13
P.....	POLIMERIZADORA	18
I.....	POLIMERIZZATRICE.....	22
D.....	POLYMERISATIONSGERÄT	27

MESTRA®



TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60


Tel. +34 944530388 - Fax +34 944711725

mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

A	Tapa / Cover / <i>Couvercle</i> / Tampa / Coperchio / <i>Deckel</i>
B	Cuadro sinóptico / Synoptic display / <i>Tableau Synoptique</i> / Quadro sinóptico / Quadro sinottico / <i>Übersichtstafel</i>
C	Teclado / Keyboard / <i>Clavier</i> / Teclado / Tastiera / <i>Tastatur</i>
D	Interruptor general / Main Switch / <i>Interrupteur Général</i> / Interruptor geral / Interruttore generale / <i>Hauptschalter</i>
E	Grifo de desagüe / Drainage tap / <i>Vanne d'évacuation d'eau</i> / Torneira de despejo / Rubinetto di scarico / <i>Ablasshahn</i>
F	Disyuntor / Circuit breaker / <i>Disjoncteur</i> / Disjuntor / Disgiuntore / <i>Ausschalter</i>
G	Contacto / Contacto / <i>Contacteur</i> / Contacto / Contattore / <i>Schütz</i>

1. Display de temperatura / **Temperature display** / *Affichage de la température* / Display de temperatura / **Display temperatura** / *Temperaturdisplay*
2. Display de tiempo / **Time display** / *Affichage du temps* / Display de tempo / **Display tempo** / *Zeitdisplay*
3. Piloto de control resistencia / **Element control lamp** / *Lampe témoin de la résistance* / Indicador luminoso de control resistència / **Spia controllo resistenza** / *Widerstandskontrollleuchte*
4. Piloto de tiempo cumplido / **Time-up lamp** / *Lampe témoin de fin de programme* / Indicador luminoso de tempo realizado / **Spia tempo trascorso** / *Zeitablaufleuchte*
5. Tecla de macha-paro / **Start-stop key** / *Digit Marche-Arrêt* / Tecla de funcionamento-paro / **Tasto avvio-arresto** / *Start-Stopp-Taste*
6. Piloto indicador de marcha / **On lamp** / *Lampe témoin de fonctionnement* / Indicador luminoso indicador de funcionamiento / **Spia di funzionamento** / *Start-Anzeigeleuchte*
7. Modo visualización de displays / **Display mode** / *Digit de changement de mode de visualitation de l'affichage* / Modo visualização dos displays / **Modalità visualizzazione display** / *Displayanzeigemodus*
8. Tecla de decrementar / **Down key** / *Digit de programmation (moins)* / Tecla de decrementar / **Tasto diminuire** / *"Verringern"-Taste*
9. Tecla de incrementar / **Up key** / *Digit de programmation (plus)* / Tecla de incrementar / **Tasto incrementare** / *"Erhöhen"-Taste*
10. Tecla de función / **Function key** / *Digit de fonction* / Tecla de função / **Tasto di funzione** / *Funktionstaste*
11. Pilotos indicadores de paso / **Step indicator lamps** / *Témoins lumineux indiquant l'étape du programme* / Indicador luminosos indicadores de passo / **Spie passo** / *Schrittanzeigeleuchten*
12. Pilotos indicadores de temperatura / **Temperature lamps** / *Témoins lumineux indiquant la température* / Indicador luminosos indicadores de temperatura / **Spie temperatura** / *Temperaturanzeigeleuchten*

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The R-080400 / R-080402 is a practical washer and polymeriser with a control panel on which sequences of polymerising times and temperatures can be programmed as required by the user. The polymeriser operates to the following 5 steps.

- Waiting time A0: a time in which the water in the tank is at ambient temperature.
- 1st ramp time A1: a time when the temperature of the water in the tank rises steadily to a final temperature set by the user.
- 1st plateau time A2: a time when the temperature of the water in the tank remains steady at the level reached at the end of the previous step. At the end of this time a warning buzzer sounds.
- 2nd ramp time A3: a time when the temperature of the water in the tank rises steadily again, to a final level set by the user.
- 2nd plateau time A4: a time when the temperature of the water in the tank remains steady at the level reached at the end of the previous step.

After these steps the unit disconnects automatically, a warning buzzer sounds and lamp (4) lights up. The time of each step can be set between 0 and 99 hours 59 minutes. The end temperatures at each plateau can be regulated by the user between 0 °C and boiling.

To enable the user to see the actual program status at all times, the polymeriser has a synoptic display. LEDs (11) show which step is ongoing and LEDs (12) also show the temperatures reached.

The unit also has two digital displays: the top display (1) shows the present temperature of the water in the tank, and the bottom display (2) shows the time remaining to the end of processing.

If W/T (9) is pressed, display (1) shows the theoretical emperature setting expected by the microprocessor, and display (2) the time remaining to the end of the step ongoing.

INSTALLATION

- Carefully unpack the polymeriser, removing the packing material inside the tank. Inside the tank you will also find a muffle tray.
- Set the unit up on a flat, horizontal, firm surface away from sources of heat or vibration. Ensure that there is a gap of at

least 15 cm between the ventilation grilles and the nearest wall or other object.

- Connect the unit to a 230 V, 50/60 Hz earthed power socket protected by a suitable thermal magnetic cut-out. Consult your installer.

PROGRAMMING

It is extremely easy to program a heating sequence for the water in the tank. Proceed as follows:

1. Press function key (10) to enter programming mode. Check that the top display reads "A0", thus indicating that this step can now be programmed. The LED for A0 (waiting time) on the synoptic display should also light up. A setting of 00.00 means that the polymeriser will start up immediately, with no waiting time.
2. Use the Up/Down keys (8 & 9) to bring up the desired first step time on the lower display (2).
3. Press function key (10). Check that the top display reads "A1", thus indicating that this step can now be programmed. The LED for A1 (1st ramp time) on the synoptic display should also light up.
4. Use the Up/Down keys (8 & 9) to bring up the desired second step time on the lower display (2).
5. Press the function key (10) again. This time the display will read "U1" to indicate that the end temperature for the first ramp can be programmed.
6. Proceed in this way until the times and temperatures for all steps have been programmed.
7. The sequence in which the steps appear on the display when the unit is in programming mode is as follows:
 - A0 - waiting time
 - A1 - 1st ramp time
 - U1 - 1st plateau temperature
 - A2 - 1st plateau time
 - A3 - 2nd ramp time
 - U2 - 2nd plateau temperature
 - A4 - 2nd plateau time
8. When programming is complete, or at any time during programming, programming mode can be exited by pressing W/T (7), which acts as an escape key. The program is memorised automatically.

9. After exiting programming mode, hold down Start/Stop (5) for approximately one second to run the polymerisation program. The program can also be stopped at any time by holding down Start/Stop (5) for around 5 seconds.

Should you wish the program to end at step A2 and omit the second ramp and plateau, merely allocate a time of 00.00 for steps A3 and A4 and the same temperature for U2 as for U1.

The temperature may not exceed 110°C at the plateaux U1 and U2. If you try to program a higher temperature, display (2) will show "A" at step U1. If you accept this, the water in the tank will be held at boiling point throughout the plateau time.

Water may boil at slightly above or below 100 °C, depending on altitude above sea level, water composition and atmospheric pressure. This means that if the boiling point programmed into the unit is 100 °C the polymerisation process may not proceed to the next step because the set temperature may never be reached.

Likewise, the polymeriser may go on to the next step when the water has not boiled properly. This problem can be solved by selecting "A" as the set temperature. The polymeriser will then automatically detect when the water boils, and hold it at that temperature throughout the step in question.

Another possibility is to detect the exact boiling point in a laboratory. The polymeriser is programmed with the following functional parameters: A0 = 00.00, A1 = 00.00, U1 = 110 °C. The temperature of the water in the tank will rise slowly to a level at which it holds steady. This is the laboratory boiling point, and this temperature should be programmed whenever you wish to polymerise at boiling point.

To prevent excessive moisture due to steam given off by the washer & polymeriser, TALLERES MESTRAITUA, S.L recommends you do not program the unit for more than 96 °C.

OPERATION

- Place the muffle tray in the bottom of the tank.
- Fill the tank with water up around 10 cm below the rim. There should always be at least 20 cm of water in the tank, to prevent damage to the heating element.
- Switch on main switch D. Some of the controls on the panel will flash, indicating that the polymeriser is ready for operation.
- Place the muffles on the muffle tray, using long, sturdy gauntlets to protect against splashes, a face-mask and muffle tongs.

TALLERES MESTRAITUA, S.L. recommends MESTRA tongs. Close over the cover.

- Program the polymeriser as per the instructions above. Once programming is completed, hold down Start/Stop for 5 seconds. The polymeriser will start up and lamp (6) will come on to indicate that the program is running as required. Lamp (3) comes on whenever the heating element is working.
- In the normal course of the program the unit can be stopped at any time by holding down Start/Stop (M/P) (5) for approximately 5 seconds. You can also enter programming mode and modify the parameters of the program which is running.
- If there is a power cut, the polymeriser retains in its memory the parameters of the step ongoing, and when power is restored the program will continue from the point at which it stopped.

If the main switch (D) is turned off before the program is completed, the unit takes this as a power cut and will therefore use maximum power when it is re-connected to attempt to retrieve the water temperature before disconnection.

Once the polymerisation program is completed, the R-080400 / R-080402 disconnects automatically and emits a warning buzz. This buzz can be cancelled by pressing W/T (7). Lamp (4) on the control panel also lights up. When removing the muffles, open the cover carefully, as steam pouring out from the tank could cause burns. Use long gauntlets and a face-mask, and remove muffles using muffle tongs. Close over the cover.

Disconnect the unit via main switch D. The fan remains on until switch D is disconnected to prevent the control module from overheating.

EXAMPLE OF PROGRAMMING

Assume that the polymeriser is to be programmed for the following sequence of steps:

1. Waiting time: zero
2. 1st ramp time: the water in the tank should reach boiling point within 1 hour.
3. 1st plateau time: the water in the tank should remain on the boil for 5 minutes. After this time the unit will buzz.
4. 2nd ramp time: the water in the tank should then drop to 80 °C in 1 hour.
5. 2nd plateau time: the water should remain at 80 °C for 30 minutes, after which the polymerisation program is completed and the unit buzzes.

The sequence of steps to be programmed is as follows:

1. Press (10) to enter programming mode. "A0" will begin to flash.
2. Press Down (8) to set a waiting time of zero.
3. Press (10) to go on to the next step. The A1 lamp will flash.
4. Press Up (9) to select a 1st ramp time of 1 hour.
5. Press (10) to go on to the next step. "U1" will flash.
6. Press Up (9) to select "A" (boiling point) as the 1st plateau temperature.
7. Press (10) to go on to the next step. "A2" will light up.
8. Press Up (9) to select a 1st plateau time of 5 minutes.
9. Press (10) to go on to the next step. "A3" will light up.
10. Press Up (9) to set a 2nd ramp time of 1 hour.
11. Press (10) to go on to the next step. "U2" will light up.
12. Press Down (8) to set a 2nd plateau temperature of 80 °C
13. Press (10) to go on to the next step. "A4" will light up.
14. Press Up (9) to select a 2nd plateau time of 30 minutes.
15. Press (10) (Function) or (7) (W/T) to memorise the program drawn up.

The polymeriser will start up when the Start/Stop button (5) is pressed for approximately five seconds, and will automatically run the program.

MAINTENANCE AND CLEANING

To obtain long, trouble-free service from the unit, follow these instructions:

- The R-080400 / R-080402 has a circuit breaker that protects the machine against electric overloads. Once the problem is solved, press the circuit breaker located on the side of the machine (F) to get it into work again.
- Message “**bro**” means that the thermal probe is broken and must be replaced.
- Message “**HHHH**” means that the triac is broken and must be replaced.
- The R-080400 / R-080402 has a circuit breaker that protects the machine against accidental overheating as it detects excessive temperature and cuts off machine operation. Once the machine returns to normal temperature, the security contactor has to be reset. To access this, proceed as follows:

1. Take off the lower panel by removing the fastening screws.
 2. Unscrew and extract the front of the machine holding the control unit.
 3. Press the contactor (G).
 4. Follow these instructions in reverse order to reassemble the unit.
 5. Carry out a test to check operation of the polymeriser.
- Clean the inside of the tank approximately every 20 hours, by heating the water to 75 °C, then emptying the tank through the drainage tap at the back.
 - Unplug the polymeriser and clean the inside of the tank with a fine scourer moistened in soapy water, taking care not to damage the heating element or temperature probe. Special products for cleaning polymerisers may also be used. Rinse the inside of the tank with plenty of water and then close the drainage tap.
 - At intervals which will depend on how hard the water is, clean away any limescale which may have been deposited on the shower nozzles or electric heating element, using a scale remover. Anti-limescale products may also be used.
 - To clean the outside housing of the unit use a cloth moistened in soapy water. Do not use solvents or other inflammable products.
 - To get the best results from the polymeriser additives available on the market to dissolve wax may be added to the water. The manufacturers' recommendations should always be followed in this case.
 - After using the polymeriser, remove wax, foam, remains of resin and other particles which may be left floating on the water.
 - For optimum performance, TALLERES MESTRAITUA, S.L. recommends the use of MESTRA muffles and flanges.
 - Clean the unit's control panel regularly with a damp cloth. Always disconnect the machine before doing this.

PRECAUTIONS

- ◆ Before connecting the unit ensure that the power supply is 230 V, 50/60 Hz and earthed. It should also be protected by a thermal magnetic cut-out. Consult your installer concerning these matters.
- ◆ Do not allow children or unskilled personnel to handle the unit.
- ◆ Use long gauntlets and muffle tongs whenever you put muffles into the tank or take them out.

- ◆ Open the cover carefully: steam pouring out from the tank could cause burns.
- ◆ Always keep the unit's cover closed while polymerisation programs are running.
- ◆ Never place more than 9 muffles on the polymerising tray.
- ◆ Always keep at least 20 cm of water in the bottom of the tank, as the heating elements may be damaged if the level is lower.
- ◆ Always leave a gap of at least 15 cm between the unit's ventilation grille and the nearest wall or other object.
- ◆ Before plugging in the unit, fill the tank with water to 10 cm below the cover.
- ◆ TALLERES MESTRAITUA, S.L. recommends polymerising at around 95 °C. At higher temperatures too much steam may be generated.

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