
## R-080406 (M18) / R-080404 (M9)

ES .....	POLIMERIZADORA CONCEPT .....	4
EN .....	POLYMERISER CONCEPT .....	8
FR .....	POLYMÉRISEUSE CONCEPT .....	11

**MESTRA®**




TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60


Tel. +34 944530388 - Fax +34 944711725

mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

A	Tapa / <b>Lid</b>
B	Panel de control / <b>Control pannel</b>
C	Interrupitor general / <b>Main switch</b>
D	Grifo de desagüe / <b>Drain tap</b>
E	Disyuntor / <b>Circuit breaker</b>
F	Contactador / <b>Contactador</b>

<b>°F</b>	Modo °F <b>Mode °F</b> <i>Mode °F</i>	<b>^</b>	Subir <b>Up</b> <i>Monter</i>
	Alarma <b>Alarm</b> <b>Alarme</b>	<b>v</b>	Bajar <b>Down</b> <i>Descendre</i>
	Programa directo <b>Direct program</b> <i>Programme direct</i>	<b>P</b>	Memoria <b>Memory</b> <i>Mémoire</i>
<b>S</b>	Calentamiento <b>Heating</b> <i>Chauffe</i>		Cambio de display <b>Display change</b> <i>Changement de display</i>
	Arranque retardado <b>Delayed start</b> <i>Demurrage différé</i>		Marcha-paro <b>Start-stop</b> <i>Marche-Arrêt</i>

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R-080404 / R-080406 is a practical polymeriser with a control panel on which sequences of polymerising times and temperatures can be programmed as required by the user.

## INSTALLATION

- Carefully unpack the polymeriser, removing the packing material inside the tank. Inside the tank you will also find a muffle tray.
- Set the unit up on a flat, horizontal, firm surface away from sources of heat or vibration. Ensure that there is a gap of at least 15 cm between the ventilation grilles and the nearest wall or other object.
- Connect the unit to a 230 V, 50/60 Hz earthed power socket protected by a suitable thermal magnetic cut-out. Consult your installer.

## ALARM

The unit allows to activate an acoustic alarm to warn when the program is finished. To activate or deactivate the alarm, proceed as follows:

- Briefly press the key . Then keep pressing it for some seconds. Pilot will blink.
- Use the keys to select “YES” or “NO” on the upper display in order to activate or deactivate the alarm.
- Briefly press the key to shift to the lower display.
- Select the type of alarm (1 or 2) by pressing the keys.
  - The alarm will go on only once.
  - The alarm will go on every 2 minutes.
- Briefly press key to confirm. Then, keep pressing it for some seconds to finish configuration.

## DELAYED START

The unit allows to program a delayed start. In order to choose a delayed start proceed as follows:

- Briefly press the key . Then keep pressing it for some seconds. Pilot will blink.
- Briefly press the key . Pilot will blink.
- Use the keys to select, on the upper display: “YES” to program a delayed start. “NO” to erase delay.
- Briefly press key to confirm.
- Use the keys to select, on the lower display, the number of hours you want the unit keeps waiting before start.

## MODE °C/°F

The unit allows to see the temperature in grades Celsius or Fahrenheit. In order to choose one mode or the other, proceed as follows:

- Briefly press the key . Then keep pressing it for some seconds. Pilot will blink.
- Briefly press the key . Pilot will blink.
- Briefly press the key . Pilot °F will blink.
- Use the keys to select, on the upper display: “YES” for °F. “NO” for °C.
- Briefly press key to confirm. Then, keep pressing it for some seconds to finish configuration.

## PROGRAMMING

In order to program a heating sequence and then save it as one of the 4 memories which can be stored, proceed as follows:

- Briefly press key .
- Simultaneously keep pressing keys . Pilot for step 1 will light.
- Use keys to select the time (minutes) to reach the first set temperatura (C1).
- Briefly press key to confirm. Pilot for step 2 will light.
- Use keys to select the set temperatura C1 that you want to reach in step 2.
- Briefly press key to confirm.
- Use keys to select the time (minutes) you want to keep temperature C1.
- Briefly press key to confirm. Pilot for step 3 will light.
- Proceed as previously to program steps 3 and 4.

If, by any reason, you want to omit one of the steps, you must simply select the value “- - -” for the time. This means that program will finish in this step.

Depending on altitude above sea level, water composition and atmospheric pressure, water may boil at slightly above or below 100 °C.

The electronic controls will automatically detect when the water boils, and hold it at that temperature throughout the step in question.

Another possibility is to detect the exact boiling point in your laboratory: set a temperature of 110 °C. The temperature of the water in the tank will rise slowly to a level at which it holds steady. This is the laboratory boiling point, and this temperature should be

programmed whenever you wish to polymerise at boiling point.

Also, to prevent excessive moisture due to steam given off by the polymeriser, TALLERES MESTRAITUA, S.L. recommends you do not program the unit for more than 95 °C. This temperature, kept for the necessary time, is enough for polymerisation.

## STORING PROGRAMS IN THE MEMORY

The electronic controls of the unit allow you to store up to 5 different memories. To store a program, proceed as follows:

- First of all, select the program you want to store, as explained in previous sections.
- Once you have selected the values you want to store, select the program (*P1*, *P2*, *P3*, *P4* or *Pr*) in which you want to store the dosage by pressing the **P** key as many times as required.
- Finally, keep the **P** key pressed for a few seconds. The display will flicker once to indicate that the values have been stored in the memory.

Memory *Pr* is reserved for the direct heating program (look at the next section).

## DIRECT HEATING PROGRAM

A direct heating program has only one heating time, only one set temperature and only one keeping time. Everytime you turn the unit on, the direct heating program will be loaded by default. This situation is indicated by pilot .

In order to modify the values of this direct heating program, simply program a cycle with only one rise and save it in memory *Pr* (look at previous section). These new values will be loaded everytime you turn the unit on. If, by error, you try to save a program with several rises in memory *Pr*, nothing will be saved.



## USING THE PROGRAMS STORED IN THE MEMORY

- To use a program stored in the memory, press the **P** key as many times as required to select the program you want (*P1*, *P2*, *P3*, *P4* or *Pr*).
- Press the key  to display the values stored.
- Press the **START/STOP**  key to start.

## OPERATION

- Place the muffle tray in the bottom of the tank.
- Fill the tank with water up around 10 cm below the rim. There should always be at least 20 cm of water in the tank, to prevent damage to the heating element.
- Place the muffles on the muffle tray, using long, sturdy gauntlets to protect against splashes, a face-

mask and muffle tongs. TALLERES MESTRAITUA, S.L. recommends MESTRA tongs. Close over the cover.

- Switch on main switch (C). Some of the controls on the panel will flash, indicating that the polymeriser is ready for operation.
- Program the unit as indicated above to ensure that the programs stored meet your own requirements. These programs can be used and modified as required in future.
- Select the program you want to run (*P1*, *P2*, *P3*, *P4* or *Pr*). If you do not select any program, then the direct heating program will run (look at the previous section).
- Press the **START/STOP**  button to start the program.
- In order to stop the running program at any time, keep pressing the **START/STOP**  button.
- Once the polymerisation program is completed, the R-080404 / R-080406 disconnects automatically. If alarm is activated, it will sound at the end of the program (look at the *Alarm* section). To remove the muffles, open the cover carefully, as steam pouring out from the tank could cause burns. Use long gauntlets and a face-mask, and remove muffles using muffle tongs. Close over the cover.
- If there is a power cut, the unit retains in its memory the parameters of the step in progress, and when power is restored the program will continue from the point at which it stopped. If the switch is turned off before the program is completed, the furnace takes this as a power cut.
- To run another program, first you must stop the unit by keep pressing the **START/STOP**  button.
- Turn the main switch off to disconnect the unit from the power supply.
- Disconnect the unit via main switch (C). The fan remains on until switch is disconnected to prevent the control module from overheating.

## DISPLAY MESSAGES

The unit can show some messages in display, to make the use easier:

- dEL** It means that the delayed start time is counting down. Lower display shows the remaining hours for the start.
- tH1** It means that the thermal probe is damaged or unplugged. Replace the thermal probe.
- tri** It means that the triac is damaged. Turn the main switch off and contact immediately to technical service.

## MAINTENANCE AND CLEANING

To obtain long, trouble-free service from the unit, follow these instructions:

- The R-080404 / R-080406 has a circuit breaker that protects the machine against electric overloads. Once the problem is solved, press the circuit breaker located on the side of the machine (E) to get it into work again.
- Message “**tH1**” means that the thermal probe is broken and must be replaced.
- Message “**tri**” means that the triac is broken and must be replaced.
- The R-080404 / R-080406 has a circuit breaker that protects the machine against accidental overheating as it detects excessive temperature and cuts off machine operation. Once the machine returns to normal temperature, the security contactor has to be reset. To access this, proceed as follows:
  1. Take off the lower panel by removing the fastening screws.
  2. Unscrew and extract the front of the machine holding the control unit.
  3. Press the contactor (F).
  4. Follow these instructions in reverse order to reassemble the unit.
  5. Carry out a test to check operation of the polymeriser.
- Clean the inside of the tank approximately every 20 hours, by heating the water to 75 °C, then emptying the tank through the drainage tap at the back.
- Unplug the polymeriser and clean the inside of the tank with a fine scourer moistened in soapy water, taking care not to damage the heating element or temperature probe. Special products for cleaning polymerisers may also be used. Rinse the inside of the tank with plenty of water and then close the drainage tap.
- At intervals which will depend on how hard the water is, clean away any limescale which may have been deposited on the shower nozzles or electric heating element, using a scale remover. Anti-limescale products may also be used.
- To clean the outside housing of the unit use a cloth moistened in soapy water. Do not use solvents or other inflammable products.

- To get the best results from the polymeriser additives available on the market to dissolve wax may be added to the water. The manufacturers' recommendations should always be followed in this case.
- After using the polymeriser, remove wax, foam, remains of resin and other particles which may be left floating on the water.
- For optimum performance, TALLERES MESTRAITUA, S.L. recommends the use of MESTRA muffles and flanges.
- Clean the unit's control panel regularly with a damp cloth. Always disconnect the machine before doing this.

## PRECAUTIONS

- Before connecting the unit ensure that the power supply is 230 V, 50/60 Hz and earthed. It should also be protected by a thermal magnetic cut-out. Consult your installer concerning these matters.
- Do not allow children or unskilled personnel to handle the unit.
- Use long gauntlets and muffle tongs whenever you put muffles into the tank or take them out.
- Open the cover carefully: steam pouring out from the tank could cause burns.
- Always keep the unit's cover closed while polymerisation programs are running.
- Never place more than 9 muffles on the polymerising tray.
- Always keep at least 20 cm of water in the bottom of the tank, as the heating elements may be damaged if the level is lower.
- Always leave a gap of at least 15 cm between the unit's ventilation grille and the nearest wall or other object.
- Before plugging in the unit, fill the tank with water to 10 cm below the cover.
- TALLERES MESTRAITUA, S.L. recommends polymerising at around 95 °C. At higher temperatures too much steam may be generated.

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