
R-080490

E.....	PULIDOR ELECTROLÍTICO.....	4
GB.....	ELECTROLYTIC POLISHER.....	6
F.....	POLISSEUR ELECTROLYTIQUE.....	8
P.....	POLIDOR ELECTROLÍTICO.....	10
I.....	PULITRICE ELETTROLITICA.....	12

MESTRA®


TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tfno. + 34 94 453 03 88 Fax + 34 94 471 17 25


E-mail: mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

1. Tapa / *Couvercle* / **Cover** / *Tampa* / Coperchio
2. Cuba de PVC / *Cuve en PVC* / **PVC tank** / *Cuba de PVC* / Vaschetta di PVC
3. Cátodo / *Cathode* / **Cathode** / *Cátodo* / Catodo
4. Ánodo (unir a la pieza) / *Anode* / **Ánodo** / *Ânodo (unir à peça)* / Anodo (unire al pezzo)
5. Amperímetro / *Ampermetre* / **Ammeter** / *Amperímetro* / Amperometro
6. Luz de puesta en marcha / *Temoin lumineux de fonctionnement* / **Start-up light** / *Luz de arranque* / Luce di avvio
7. Fusible / *Fusible* / **Fuse** / *Fusible* / Fusibile
8. Potenciómetro de ajuste / *Potentiomètre* / **Regulating** / *Potenciómetro de regulação* / Potenziometro di regolazione
9. Reloj temporizador / *Minuterie* / **Timer** / *Relógio temporizador* / Timer
10. Borne positivo / *Borne (+)* / **Positive terminal** / *Borne positivo* / Morsetto positivo
11. Borne negativo / *Borne (-)* / **Negative terminal** / *Borne negativo* / Morsetto negativo
12. Disyuntor 10 A / *Disjoncteur 10 A* / **Circuit breaker** / *Disjuntor 10 A* / Disgiuntore 10 A

1

A: Ánodo (+), superficie del metal a pulir.
Ánode (+), surface of the metal.
Anode (+), surface du métal à polir.
 Ânodo (+), superfície do metal a polir.

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The R-080490 is a practical electrolytic polisher especially designed for use in jewellery work and dental mechanics. So that it may work properly and perform polishing of high quality. We recommend that you carefully read the following instructions.

OPERATION

Any metal studied under the microscope exhibits a more or less wrinkled surface, with projections of every sort. If we put the same metal into an electrolytic bath (fig. 1) and connect it to a component with positive electric voltage (anode), there is a curious effect: the electrolyte withdraws, giving rise to charged molecules, called ions, that tear atoms off the surface of the metal to combine with them. The resulting particles have an electric charge, for which reason, when they are submitted to the difference in potential between cathode and anode (fig. 5), there is a smooth displacement toward the negative pole (cathode). The metal particles are deposited on the cathode (fig. 3), and there is again a charged ion that starts the cycle off anew. The process comes to an end when the electric flow is cut.

The parts of the metal surface on which there appear the most pronounced projections are those that have a large surface contact with the electrolyte, hence those that lose most metal atoms during polishing. In consequence, when polishing is over the surface of the metal has much less of an escarpment appearance (fig. 4), with the rough edges smoothed down, so that light can be reflected more homogeneously. The metal is therefore brighter. The excess of metal removed from the anode (+) ends up deposited on the surface of the cathode.

USE OF THE UNIT

1. Position the apparatus on a surface that is flat, firm and fully horizontal. There are four rubber feet on the base to ensure that it will stay put.
2. Place the PVC tub (2) next to the unit. Fill the tub with about 1.5 litres of electrolyte. Be careful to ensure that the electrolyte is of high quality, for in great measure the efficiency of the polishing will depend on this factor.
3. Connect the wiring between tub and unit. In doing so you should take into account the colour coding. BLACK corresponds to the cathode and should be connected to the annular electrode at the back of the tub (3). RED corresponds to the anode and should be connected to the copper rod on top of the tub (4).
4. To the tub anode, connect the metal parts that you wish to polish. For this purpose, the unit is supplied with a coil of silver-plated copper wire. Cut a few small sections of wire and use them to make hooks that will allow you to hang, from the tub (copper red on the top), the metal pieces that are to be polished. Take care to ensure that the electrical contact with this electrode is sound and that the length of the hook allows, in addition, the piece to be fully submerged in the bath. Then close the tub cover (1).
5. Insert the plug into a 230V, 50/60 Hz outlet with earth.
6. Start up the unit via the main switch on the regulation timer (8). Check that the little pilot light (6) goes on. Adjust the current of the electrolytic bath by means of the potentiometer. The intensity of the instantaneous current is shown on the ammeter (5). To proceed with polishing, it is advisable to use 0,5 amperes per cm of metal surface. Since unit current should not exceed 10 A, the surface for polishing should never exceed 20 cm². For orientation, the surface of a skeletal is usually between 6 and 10 cm². Ibis means a current between 3 and 5 amperes per skeletal, hence a maximum of 2-3 skeletal per bath.
7. The unit has a practical timer (9). The period over which the bath will function can conveniently be set anywhere up to 30 minutes. If polishing is done according to the specifications regarding current in the foregoing point, time should usually be between 10 and 15 minutes. This aspect of the matter will in any case depend on the current chosen.
8. It is probable that during polishing the electrolyte in the tub will heat up. This is perfectly normal. Also it is possible that by accident those pieces for polishing that are positioned on the anode hooks touch the annular cathode of the tub. When this occurs there is a short circuit and the circuit breaker (12) trips automatically. In such circumstances, remove the pieces causing the short circuit and wait three minutes before resetting the circuit breaker.

9. Another cause of circuit breaker disconnection may be excess of current in polishing. If there are a lot of pieces, current may exceed maximum values. When this happens, remove the excess pieces and wait three minutes before resetting the circuit breaker.
10. At the end of polishing time, the unit will disconnected automatically. Remove the polished pieces and wash them in abundant water. If by mistake you have connected the wrong polarities, you will find that the piece is coated with a fine layer of copper. Reconnect cathode and anode correctly and repeat the process in order to remove the copper layer.

MAINTENANCE AND CLEANING

To ensure that the unit functions at its best and has a long life, we recommend the following:

- Always use electrolytes of recognized quality. Ibis is especially important if you want optimum polishing.
- Regularly renew the electrolyte according to the maker's instruction.
- So that sound electrical contact is ensured, frequently renew the wire hooks that hold the pieces. Talleres Mestraitua, S.L. recommends that you use the wire **MESTRA**.
- The anode is the copper rod on top of the tub, from which the pieces hang that are to be polished to remove oxide. Clean it regularly. Use fine paper.
- Regularly empty the tub and remove its electrodes in order to clean thoroughly. Use a scourer and soap.
- Do not let the electrolyte come in contact with the control panel. To clean the surface of the panel, use a cloth dampened in soapy water.
- Do not overload the unit. It is recommended that operating current be kept below 7.5 A. In addition, the surface area of a piece for polishing should not exceed 15 cm² (2-3 skeletal).
- Avoid short circuits between the electrodes of the unit. If the circuit breaker trips, wait at least three minutes before resetting it.

PRECAUTIONS

- ◆ Make sure that manis voltage is correct (230 V, 50/60 Hz, and earthed).
- ◆ Do not use the unit in excessively humid places.
- ◆ See to it that the unit is not used by children or by personnel not specialized.

- ◆ Before cleaning it, make sure the unit is unplugged.
- ◆ Do not submerge the unit in any liquid for cleaning purposes.
- ◆ Do not put any liquid into the tub other than an electrolyte of high quality.
- ◆ While the unit is operating, the electrolyte will probably heat up. Certain fumes will result. Do not inhale them.
- ◆ Carefully read the instructions of the electrolyte manufacturer. Replace the electrolyte as often as specified. Handle it with care to avoid spillage. In addition, avoid prolonged skin contact.
- ◆ To dispose of the electrolyte when it has become degenerate, follow the instructions of the manufacturer in relation to respect for the environment
- ◆ Do not overload the unit. Ensure that operating current does not exceed 7,5 amps.
- ◆ Keep the top of the tub closed as long as possible, for in this way you will prevent uncontrolled spillage of electrolyte.

TECHNICAL CHARACTERISTICS

Voltage:	230 V, 50/60 Hz
Secondary voltage:	16 V
Maximum current:	10 A
Circuit breaker:	10 A
Timer:	0-30 minutes
Long x wide x high:	245 x 245 x 140 mm
Weight:	5.3 kg
Tank capacity:	2 litres
Tank height x diam.:	160 x 145 mm
Tank weight:	1.2 kg

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