**R-080511
ANETO**

E.....	MEZCLADORA DE GELATINA	4
GB.....	GEL MIXER.....	8
F.....	MACHINE A GELATINE.....	12

MESTRA®

TALLERES MESTRAITUA S.L.

Txori-Erri Etorbidea, 60

Tfno. + 34 944530388 - Fax + 34 944711725

E-mail: mestra@mestra.es - www.mestra.es

48150 SONDIKA - BILBAO - ESPAÑA

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The R-080511 is a practical gel mixer specially designed for use within the mechanical dentistry sector. It enables the gel placed in the container to a heating and stirring lea defined by the user which consists of three phases:

1. A prior waiting period, during which the unit remains on standby awaiting the heating phase to commence.
2. A heating period, in wich the gel temperature rises to a value (C1) set beforehand by the user, between 65 and 95 °C. During this phase heating temperature is higher than the subsequent pouring temperature. The idea is to achieve a homogeneous gel mixture and without any clots. Also, this phase the mixer blades are set in motion in order to further favour the correct mixing of the gel.
3. Once the heating phase the gel achieves an optimum texture, but its temperature is too high for pouring. A cooling period then commences forced by a fan which ends when the gel temperature reaches a value (C2) preset by the user (30-65 °C). Once the third phase has been completed, the gel is then ready for pouring

INSTALLING THE UNIT

- Carefully unpack the gel mixer.
- Choose a flat surface to place it on, far from any heat or vibration sources. We also recommend choosing a location which is sufficiently ventilated.
- Be sure to leave a space between the surface of the table and the legs top, given that the cooling air from the container is collected. Also, leave 15 cm between the back cooling grid and the wall.
- Conect the unit to a 230 V, 50/60 Hz electric outlet with an earth connection.

TEMPETARUTE SETTING

1. Connect the main switch of the unit (F).
2. Upper and lower displays will show the C1 and C2 temperatues respectively.
3. You can modify both of them by pressing the **▲▼** keys.
4. To shift from a display to the other press the **☐** key. The doth at the rightside indicates the active display.
5. Both temperaturas can be modified even when the unit is working.

ALARM

The R-080511 Gel Mixer allows you to activate an acoustic alarm to warn when the gel is ready to pour. When the unit finishes the cycle, the display will read "End". If the alarm is on (🔔 pilot light on), an acoustic signal will go off as well.

To activate or deactivate the alarm:

- Press the **☐** key for a few seconds to entering the configuration mode.
- You will see that the 🔔 (alarm) pilot light is flickering.
- Select "YES" or "NO" on the upper display in order to activate or deactivate the alarm by pressing the **▲▼** keys.
- Press shortly the **☐** key to shift to the lower display.
- Select the type of alarm (1 or 2) by pressing the **▲▼** keys.
 - 1 The alarm will go on only once.
 - 2 The alarm will go on every 2 minutes.
- Press the **☐** key for a few seconds to escape from the configuration mode.

DELAYED STARTUP

The R-080511 Gel Mixer allows you to program a delay startup:

- Press the **☐** key for a few seconds to entering the configuration mode.
- You will see that the 🔔 (alarm) pilot light is flickering.
- Press shortly the **☐** key to see at the upper display the message "**dEL**" (delayed startup).
- By pressing the **▲▼** keys, select at the lower display the time (hours) you want the unit keeps waiting before start heating.
- Press the **☐** key for a few seconds to escape from the configuration mode.
- Now, when you press the **START/STOP** button (C) to start the unit, upper display will show the message "**dEL**" and the lower one will show the number of resting hours. After the waiting time, the unit will start up.

TEMPERATURE IN °C/°F

The R-080511 Gel mixer allows you to read the temperatura in °C or °F:


- Press the **☐** key for a few seconds to entering the configuration mode.
- You will see that the 🔔 (alarm) pilot light is flickering.

- Press shortly the  key several times to light the icon °F (Fahrenheit degrees).
- By pressing the ▲▼ keys, select at the upper display the value “YES” to work in °F or “No” to work in °C.
- Press the  key for a few seconds to escape from the configuration mode.
- Icon °F lighting indicates that temperatura is shown in °F.

OPERATION

1. Open the upper lid (A).
2. Fill the container with gel. In order to do this, cut the solid block of gel into small pieces of about 4 cm long, and place them in the container (do not cut up the gelatine inside the tank). The final level of the gel should be about 3 cm, above the upper blade. When the unit enters into operation, the level of the gel in the container will fall somewhat due to the heating resistance and blades. Then refill the container again until the final level of the gel covers the upper blade.

WATCH OUT!! Do not let solidify the jelly into the vessel: when the machine start up, the temperature sensor could take a wrong value and consequently a wrong machine perform. If this happens, full empty the vessel, cut the jelly in small pieces and start up again the machine.

3. Close the lid and programme the unit following the previously indicated instructions. Press the **STOP/START** button (C) to start the cycle. The unit will startup automatically if a delayed startup has not been programmed. During the operation, the pilots  (resistance heating) and  (motor working) will light in order to indicate to the user that the programme continues to run without any problems. The upper display will show the temperatura inside the tank. The lower display will show C1 or C2 to indicate if the unit is in the heating or cooling phase.
4. To stop the operation of the unit at any moment keep the **START/STOP** button (C) pressed.
5. In case a delayed startup is programmed, upper display will show the message “dEL” and the lower one will show the number of resting hours. After the waiting time, the unit will start up.
6. You can modify the values of the C1 and C2 temperatures at any time during the programme execution.
7. Once the heating temperature (C1) is reached, the the cooling phase will start until the pouring

- temperature (C2) is reached. The lower fan will work. When this phase finish, the display will show the message “End”.
8. When the pouring temperatura (C2) is reached, gel is ready to work. Place a container under the spout (E) and pull the knob (F) towards you. The gel will flow out.
9. Once the work with the unit has been completed, empty the unit container by pouring all of the excess gel into another container through the spout (E).
10. In order to finish the programme, keep pressing the **START/STOP** button (C) for some seconds. You will notice that the icons  and  start blinking. Then disconnect the main switch.

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During the normal running of the programme the lid of the unit cannot be opened. A safety switch will disconnect the unit if the upper lid (A) is opened, and the display will show an error message “dOO”.

STORING PROGRAMS IN THE MEMORY

The electronic controls allow you to store up to 4 different programs. To store a program, proceed as follows:

- First of all, set the two tempetarues C1 and C2 you want to store, as explained in previous sections.
- Once you have selected the values you want to store, select the program (P1, P2, P3 or P4) in which you want to store the dosage by pressing the **P** key as many times as required.
- Finally, keep the **P** key pressed for a few seconds. The display will flicker once to indicate that the values have been stored in the memory.

USING THE PROGRAMS STORED IN THE MEMORY

- To use a program stored in the memory, press the **P** key as many times as required to select the program you want (P1, P2, P3 or P4).
- Press the key  to display the values stored.
- Press the **START/STOP** button (C) to start.

POWER SUPPLY BREAKDOWN

If power supply is interrupted (or main switch is disconnected) while the unit is working, once the power is recovered, the electronic control will try to continue the program at the point it was stopped.

In order to correctly stop the unit, keep pressing the **START/STOP** button (C) for some seconds until pilots  and  blink, before disconnecting the main switch.

BLOCKED MOTOR

If the electronic control detects that the motor is blocked (by example by solid gel), the display will show the error message "cte".

The user will have to choose between two options:

1. Keep pressing the **START/STOP** button (C) to end the cycle and solve the problem
2. Briefly press the **START/STOP** button (C) to allow the unit heat the gel for 10 minutes before trying to start the motor again.

MAINTENANCE AND CLEANING

In order to guarantee optimum reliability and long working life of the unit, we recommend you bear in mind the following:

- Do not use solvents or any other inflammable product.
- The unit is provided with a thermoprobe. If these were to cause any problems, immediately contact the **MESTRA** technical service.
- The motor reducer which drives the blades of the stirrer is protected by a thermal fuse and by the microprocessor's own programming software, which disconnects the motor whenever overloads occur.
- All the units are checked before leaving the factory. During the rigorous final quality control carried out in TALLERES MESTRAITUA, S.L. several test are made using gel, in order to better simulate the true use conditions of the unit. For this reason, it is completely normal to receive your unit with some remains of gel inside the container.

PRECAUTIONS

- ◆ **Before carrying out any maintenance on the unit, switch it off and unplug from the power supply.**
- ◆ Before connecting the unit make sure the electric connection is at 230 V, 50/60 Hz with an earth connection.
- ◆ Neither children nor unqualified personnel are allowed to manipulate the unit.

- ◆ Do not open the upper lid when the unit is in operation. Under circumstances must you bridge the safety switch which protects the opening of the lid.
- ◆ Avoid the level of gel rising above 3 cm over the last blade (the closest to the upper edge of the container).
- ◆ Place the unit upon a flat surface, correctly levelled, far from any sources of vibration. Avoid the possibility of it tipping over.
- ◆ Make sure the tap is turned off before starting up the unit.
- ◆ Do not pour the gel at a temperature (C1), because it may scald you. Wait until the pouring temperatures reaches (C2).
- ◆ During operation avoid touching the lid. It may have a high temperature. In the gel pouring manoeuvre, wear gloves for protection.
- ◆ Be extremely careful during the start ups after a motor blocking or power supply breakdown. Eventual spattering or spillage of gel may occur.
- ◆ Never leave the gel in liquid state during a delayed start up. It may solidify, preventing a normal start up of the unit. Always insert the gel in pieces when programming a delayed start up.

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