

Cenit L



Technical Datasheet

Product description: Light-curing blocking-out resin.

Manufacturer: Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

Composition:

Pos	chemical designation
1	aliphatic urethanacrylate
2	tetrahydrofurfuryl-2-methacrylate

Properties:

Colour: light blue
Spectral range: 300 – 500nm
Vickers hardness: > 600 N/mm²
Flexural strength: > 80 N/mm²
E-modulus: > 2,0 kN/mm²

Dreve Dentamid GmbH
Max-Planck-Str. 31
59423 Unna / Germany
Tel.: 02303 / 8807-0
Fax: 02303 / 8807-55
info@dreve.de
www.dreve.com