



MESTRA[®]

DECANTADORA CON BOLSA (Tamaño reducido)
DECANTER WITH COLLECTION BAG (Small size)
BAC DE DÉCANTATION AVEC SACS (Petit modèle)
R-080087



R-080085-02
Bolsas de repuesto (20 unidades)
Spare bags (20 units)
Sacs de remplacement (20 unités)



MESTRA®

DECANTADORA CON BOLSA (Tamaño reducido)
DECANTER WITH COLLECTION BAG (Small size)
BAC DE DÉCANTATION AVEC SACS (Petit modèle)
R-080087

Nº	Ref.	DESCRIPTION	DESCRIPTION
1	080087-01	Cuba	Container
2	080087-02	Tapa	Lid
3	080085-08	Casquillo fijación cartabones	Squares fixation sleeve
4	080087-11	Cartabones	Squares
5	080087-05	Tubo recto entrada	Inlet straight tube
6	080087-06	Adaptador salida desagüe	Outflow exit adaptation
7	080087-07	Tubo flexible salida	Flexible exit tube
8	080087-08	Abrazadera	Hose clamp
9	080087-09	Tapón roscado	Threaded plug
10	080085-13	Cadena	Chain

Bolsas de repuesto
Spare bags
Sacs de remplacement

11	080085 - 02	Bolsas de repuesto (20 U.)	Spare part bags (20 U.)
----	-------------	----------------------------	-------------------------