


MESTRA®
 MEZCLADORA DE VACIO (KOALA)
 (KOALA) VACUUM MIXER
 MALAXEUR SOUS VIDE (KOALA)
 R-080520



MESTRA®

MEZCLADORA DE VACIO (KOALA) (KOALA) VACUUM MIXER MALAXEUR SOUS VIDE (KOALA) R-080520

Nº	Ref.	Concepto	Description	Nº	Ref.	Concepto	Description
1	080520-01	Carcasa	Housing	19	080240-31	Filtro silenciador	Silent filter
2	080520-02	Carátula control	Control panel	20	080520-19	Retén	Seal
3	080520-03	Control	Control	21	000336	Racord recto metálico 1/8 Ø6	1/8 straight Ø6 connector
4	080520-04	Vacuómetro	Vacuum meter	22	080520-21	Pipa	Pipe
5	000333	Racor recto hembra	Female straight conector	23	080520-45	Llave de vaso	Drive socket
6	080520-06	Bomba de vacío	Vacuum pump	24	080518-19	Tornillo llave de vaso	Drive socket screw
7	080520-07	Motorreductor	Gear motor	25	080520-24	Junta tórica de pipa	Pipe O-ring
8	080520-08	Distribuidor	Distributor	26	000250	Cable	Cable
9	000326	Racor codo metálico 1/4 Ø8	1/4 elbow Ø8 connector	27	080523-08	Tuerca de eje	Shaft nut
10	000312	Racor codo plástico 1/8 Ø6	1/8 elbow Ø6 connector	28		Tapa de taza	Bowl lid
11	000308	Racor codo plástico 1/8 Ø4	1/8 elbow Ø4 connector	29		Junta tórica de taza	Bowl O-ring
12	000310	Racor codo plástico 1/4 Ø6	1/4 elbow Ø6 connector	30		Aspas de taza	Bowl blades
13	080520-13	Chasis	Support	31		Vaso de taza	Bowl cup
14	080520-14	Chapa sujeta filtro	Filter holder sheet	32	080520-26	Pletina sujeción máquina	Machine support sheet
15	080520-15	Filtro	Filter	33	080520-27	Angulo sujeción máquina	Machine support angle
16	080520-16	Base con interruptor	Plug base with switch	34	080520-28	Pletina sujeción pared	Wall support sheet
17	080520-17	Condensador	Capacitor	35	000231	Fusible 4A	Fuse 4A
18	080242-23	Electroválvula	Electrovalve	36	060002	Espátula de plástico (12 U.)	Plastic spatula (12 U.)

Bowl 250 cc. R-080523	
27	080523-08
28	080523-01
29	080523-04
30	080523-02
31	080523-03

Bowl 550 cc. R-080524	
27	080523-08
28	080524-01
29	080524-04
30	080524-02
31	080524-03

Bowl 750 cc. R-080525	
27	080523-08
28	080525-01
29	080525-04
30	080525-02
31	080525-03

Bowl 1000 cc. R-080526	
27	27-080523-08
28	28-080526-01
29	29-080526-04
30	30-080526-02
31	31-080526-03

Bowl 1250 cc. R-080527	
27	080523-08
28	080527-01
29	080527-04
30	080527-02
31	080527-03

Bowl 230 cc. R-080528	
27	080523-08
28	080523-01
29	080528-04
30	080528-02
31	080523-03