


Spare parts 18080500 Version A



Reparaturzeichnung • Repair services drawings • Eclatés pour les services de réparation • Esplosi per servizi d'assistenza
 Diseños para servicios de reparaciones • Загрузка чертежей Для ремонтных мастерских

1808 xxxx MT3



21-6630 14122021 / A

Position	Product name	Item number
1	Admission hose compl. Trimmer	900020976
2	Outlet hose 41x35mm grey 1,5m long	900150449 *
3	MT3/MT3 pro Trimmer table complete	900036292
4	Screw for trimming disc MT2/TT2/MTplus/MT3/MT3 pro	900035237
5	Door locker MT3/MT3 pro	900036267
6	Spray pipe MT3/MT3 pro	900036405 *
7	Hardware for MT3/MT3 pro	918080001
8	Tools for MT3/MT3 pro	900036614



Position	Product name	Item number
9	Flange MT3/MT3 pro complete	900036296
10	Set of O-rings for flange MTplus/MT2/TT2/MT3/MT3 pro	900033450 *
12	Door gasket MT3/MT3 pro	900036511 *
13	Rubber bases	900034653
14	Door complete MT3	900036286
15	Cap on/off switch MT3	900036412 *
16	Bracket MT3	900036293
17	Holder for trimmer table MT3 / MT3 pro	918080004
101A	Driving motor MT3 230V	900036282
101B	Driving motor MT3 120V	900036321
102A	Water admission 230V MT3 / MT3pro	918080002
102B	Water supply 100/120V MT3 / MT3pro	918080003
103A	Engine condenser 16 μ F	900134609
103B	Engine condenser 50 μ F MT3/MT3 pro	900033457
104	Joint f.split ring MT3/MT3pro	900036612
105	Reed switch MT3	900036275
106A	On/Off-switch MT3 230 V/12 A	900138698
106B	On/Off-switch 120V MT3	900135736
106C	On/Off-switch MT3 230 V/10 A	900033777
107A	Power cord MT3/MT3 pro 230 V	900036260
107B	Power cord MT3/MT3 pro 120 V	900036262
107C	Power cord MT3/MT3 pro / UK	900036967
107D	Power cord MT3/MT3 pro / China	900036968
107E	Power cord MT3/MT3 pro / Brazil	900036969
107F	Power cord MT3/MT3 pro / South Africa	900036970
107G	Power cord MT3/MT3 pro / Suisse	900036971

Position	Product name	Item number
107H	Power cord MT3/MT3 pro / Argentina	900036975
107I	Power cable MT3 (230V-10A) / AUS/NZL	900033778
501	Klettfix Dia.234mm MT3/MT3pro	18031001
502	Klettfix MTplus, MT3 and MT3 pro, refill pack grain 80	18031100
503	Klettfix MTplus, MT3 and MT3 pro, refill pack grain 120	18031200
504	Marathon MT3/MT3pro Dia. 234mm	18032001
505	Infinity MT3/MT3pro Dia. 234mm	18033001

* Wear and tear parts, excluded from the guarantee.

