


Spare parts 18460000 Version A



Ersatzteilzeichnung • Spare parts drawings • Plan des pièces de rechange • Disegno delle parti di ricambio
Plano de piezas de recambio • Чертеж запасных частей

1846*000 POWER Steamer 2



21-6924 27082025_A

Position	Product name	Item number
1	Tank cap, cpl.	900021865
2	Service flap	900033741
3	Tool (socket wrench)	900033705
4	Set of seals	900021070 *
5	Solenoid valve steam	900037875 *
6	Rubber bases	900034653
7	Handpiece holder cpl.	900033666
8	Cap (service opening)	900033642
	Front part, center, PS2	900033639



Position	Product name	Item number
10	Nozzle	900037862 *
11	Complete handpiece kit	900021682
12	Admission hose compl. Trimmer	900020976
13	Cleaning brush	900033800
101	Pressure switch, steam, cpl.	900021880
102A	3 pin power cord	900033712
102B	3 pin power cord GB	900038903
102C	3 pin power cord Argentina	900038909
102D	3 pin power cord Brazil	900038905
102E	3 pin power cord South Africa	900038906
102F	3 pin power cord AUS/NZL	900038910
102G	3 pin power cord USA	900038900
102H	3 pin power cord Switzerland	900038907
102I	3 pin power cord China	900038904
103A	Circuit board (control) PS2 230 V	900033625
103B	Circuit board (control) PS2 120 V	900033624
103C	Circuit board (control) PS2 100 V	900033623
104	Circuit board (operation)	900033626
105	Front part, top, incl. key pad PS2	900033750
106	Seal for connecting piece	900033651
107	Support blocks set	900033663
108	Temperature sensor	900033709
109	Boiler thermal fuse 298°C	900030258
110	Pressure relief valve	900021927
111	Micro switch single	900033603
112	Transformer for voltage variants	900033671
113	Isolation set	900021751

Position	Product name	Item number
114	Solenoid valve water	900037874
115	Pressure switch, water cpl.	900021881
116	Pump	900021882
117	Ribbon cable	900033728
118	Check valve	900037878
119	Chassis, right side	900033637
120	Chassis, left side	900033638
121	Chassis, upper part	900033640
122	Pressure gauge incl. screw connection	900021892
123	EPDM foam hose	900033937
124	Plastic base POWER steamer cpl.	900030276
125A	Boiler POWER steamer 2 230V	900021755
125B	Boiler POWER steamer 2 120V	900021745
125C	Boiler POWER steamer 2 100V	900021748
501	POWER steamer Wall bracket	18450200

* Wear and tear parts, excluded from the guarantee.

