


Spare parts 21570000 Version A



Reparaturzeichnung • Repair services drawings • Eclatés pour les services de réparation • Esplosi per servizi d'assistenza
 Diseños para servicios de reparaciones • Загрузка чертежей Для ремонтных мастерских

2156 xxxx / 2157 xxxx Waxlectric I / II



21-9670 2311/2023 / A

Position	Product name	Item number
1	Heating handle, red	21540002 *
2	Heating handle, yellow	21540001 *
3	Cork-grip	900035948 *
4	Power supply	900033752
5	Retainer ring	900032962 *
6	Rubber bases	900034564
7A	Power cord Waxlectric Europe	900035975
7B	Power cord Waxlectric USA	900035976
7C	2-pin power cord / UK	900038652



Position	Product name	Item number
7D	2-pin power cord / China	900038653
7E	2-pin power cord / Brazil	900038654
7F	Power cord – Argentina	900038680
7G	2-pin power cord / AUS/NZL	900032773
101A	Mains circuit board Waxlectric I	900036125
101B	Mains circuit board Waxlectric II	900036126
102A	Operating foil Waxlectric I	900036127
102B	Operating foil Waxlectric II	900036128
501	Kit of holders Waxlectric	21511500
502	Probe small for WAXLECTRIC	21550101 *
503	Probe medium for WAXLECTRIC	21550102 *
504	Probe large for WAXLECTRIC	21550103 *
505	Blade narrow for WAXLECTRIC	21550104
506	Blade wide for WAXLECTRIC	21550105
507	Hollenback horizontal for WAXLECTRIC	21550107
508	Hollenback vertical for WAXLECTRIC	21550108
509	Probe mini angled for WAXLECTRIC	21550109
510	Lance shaped blade for WAXLECTRIC	21550110
511	Beaver tail blade for WAXLECTRIC	21550111
512	Wax knife large for WAXLECTRIC	21550112

* Wear and tear parts, excluded from the guarantee.