


Spare parts 22000400 Version A



Reparaturzeichnung • Repair services drawings • Eclatés pour les services de réparation • Esplosi per servizi d'assistenza
Diseños para servicios de reparaciones • Загрузка чертежей Для ремонтных мастерских

2200 0400 Mobiloskop 'S'



21-5802 19012018 / A

Position	Product name	Item number
2	Eye-piece cup for Microscope soft	922000010 *
3A	Pair of oculars, 10x 2 pcs.	22000110
3B	Pair of oculars, 5x 2 pcs.	22000105
4	Ocular holder for microscope, adjustable	922000003
5	Protective glass for Mobiloskop objective lenses	900030743 *
6	Protective cover for Mobiloskop	900037683
7	Retainer Ring 2200 2 pcs.	900030737 *
8	Fastening adapter for Mobiloskop	900021573



Position	Product name	Item number
9	Support arm with universal clamp	22000401
11	Stand with screw clamp	22000420
12	Mobiloscope LED lighting 100-240 V, 50/60 Hz	22005000

* Wear and tear parts, excluded from the guarantee.

