




Reparaturzeichnung • Repair services drawings • Eclatés pour les services de réparation • Esplosi per servizi d'assistenza
 Diseños para servicios de reparaciones • Загрузка чертежей Для ремонтных мастерских

2934 0000 / 2934 1X00 SILENT compact



21-6830 22022021 / A



Version	Position	Name	Article number	Wishlist
A	1	Exhaust air filter kit	900037585*	
A	2	Fine filter SILENT compact / Silent compact CAM	900021635*	
A	3	Bellow seal	900037473*	
A	4	Turbine SILENT compact	900037226*	
A	5A	Mains cable 230 V with low-power plug	900032750	
A	5B	Mains cable 120V with low-power plug	900035570	
A	5C	Unit connector cable 230V UK	900035938	
A	5D	Power cord Argentine with cable-mount connector, straight	900032769	
A	5E	Unit connector cable 230V CH	900036070	
A	5F	Power cord Brazil with cable-mount connector, straight	900032767	
A	5G	Power cord SA 10A with cable-mount connector, straight	900032765	
A	5H	Power cord Italy with cable-mount connector, straight	900032768	
A	5I	Unit connector cable 230V OZ lightgrey	900035840	
A	6	Kit of small parts	900037238	
A	7A	Schuko adapter / UK	900033552	
A	7B	Schuko adapter / China	900033553	
A	7C	Schuko adapter / Suisse	900033556	
A	7D	Schuko adapter / South Africa	900033555	

*) Wear and tear parts, excluded from the guarantee.

Version	Position	Name	Article number	Wishlist
A	7E	Schuko adapter / Italy	900033557	
A	8	Silent dust drawer	900037230	
A	101A	Printed board circuit SILENT compact	900037408	
A	101B	Printed board circuit SILENT compact, 120V	900033513	
A	102	Operating foil SILENT compact	900037234	
A	103	Appliance protective switch 8A	900037389	
A	104A	Toroidal core transforme 230VAC / 120VAC	900037224	
A	104B	Toroidal core transformer 220 VAC / 100 VAC	900037682	
A	105	Foam inlay complete	900021590	
A	106	Vibration motor for cleaning	929349998*	
A	501	End bushing set, 2 pieces	29210003	
A	502	90° angled connector Silent	29340007	
A	503	Suction Hose	900034240	
A	504	Suction hose, antistatic, 3 m, incl. 2 end mufflers	900034826	
A	505A	Suction hose inside diameter 38 mm, 6 m	901150823	
A	505B	Suction hose inside diameter 38 mm, 9 m	902150823	
A	506	Adapter for hose connection	900034305	
A	507	Universal suction hose adapter	900034430	

*) Wear and tear parts, excluded from the guarantee.

Replacement parts

SILENT compact, 220-240 V



Version

A

Article number

29340000



making work easy

Version	Position	Name	Article number	Wishlist
A	508	External exhaust air adapter SILENT compact	29340004	
A	509	Extractor clamp	29250000	
A	510	Glass pane with attachment	29251000	

*) Wear and tear parts, excluded from the guarantee.