




Reparaturzeichnung • Repair services drawings • Eclatés pour les services de réparation • Esplosi per servizi d'assistenza
 Diseños para servicios de reparaciones • Загрузка чертежей Для ремонтных мастерских

2936 X000 SILENT TC2



21-6838 18012019 / A



Version	Position	Name	Article number	Wishlist
A	1	SILENT TC, TC2 Fine filter Class M	900021532*	
A	2	Exhaust air filter set SILENT EC2, powerCAM EC	900037625	
A	3	Air filter for electronics complete	900034224*	
A	4	Retainer rings set SILENT EC2, powerCAM EC	900037616*	
A	5	Motor unit SILENT TC, TC2, 230 V	900021631*	
A	6	Hardware kit SILENT EC2, powerCAM EC	900038671	
A	7	Pinch valve, straight SILENT TS2/ TC2 / EC2 / V4	900021526*	
A	8	Air intake duct for SILENT TC/ TC2/ EC2/ powerCAM EC	29370001*	
A	9	Pneumatic connection set	21068	
A	10	Dust drawer SILENT EC2, powerCAM EC	900037423	
A	11	Vibration motor for cleaning SILENT EC2, powerCAM EC	900021589	
A	12	Filter nut	900037757	
A	13A	Schuko adapter / UK	900033552	
A	13B	Schuko adapter / China	900033553	
A	13C	Schuko adapter / South Africa	900033555	
A	13D	Schuko adapter / Suisse	900033556	
A	14	SILENT TC, TC2 Exhaust air labyrinth with metal grid	900021695*	
A	101A	Power cable 230V	900033425	

*) Wear and tear parts, excluded from the guarantee.

Version	Position	Name	Article number	Wishlist
A	101B	Power cord 120V	900033426	
A	101C	Power cord UK	900033543	
A	101D	Power cord Millo /pro / China	900033544	
A	101E	Power cord Millo /pro / Brazil	900033545	
A	101F	Power cord Millo /pro / South Africa	900033546	
A	101G	Power cord Millo /pro / Suisse	900033547	
A	101H	Power cord Argentine	900033560	
A	101I	Power cord 3x1,5qmm / AUS/NZL	900036979	
A	103A	SILENT TC2 Printed circuit board 230 V	900037433	
A	104	Control unit cpl. SILENT TC2	900021600	
A	105	Appliance protective switch 15A	900034668	
A	106A	SILENT TC, TC2 Toroidal core transformer	900037411	
A	107	Pneumatic unit Silent TS2	900036814	
A	501	Waste bag for TC, TC2, EC2 & powerCAM EC (25 pcs)	29340014*	
A	502	Waste bag for SILENT TC, TC2, EC2 Starter Kit	29340015*	
A	503	End bushing set, 2 pieces	29210003	
A	504A	Suction Hose	900034240	
A	504B	Suction hose, antistatic, 3 m, incl. 2 end mufflers	900034826	

*) Wear and tear parts, excluded from the guarantee.

Replacement parts

SILENT TC2, 220-240 V



Version

A

Article number

29360000



making work easy

Version	Position	Name	Article number	Wishlist
A	505A	Suction hose inside diameter 38 mm, 6 m	901150823	
A	505B	Suction hose inside diameter 38 mm, 9 m	902150823	
A	506	Adapter for hose connection	900034305	
A	507	Universal suction hose adapter	900034430	
A	508	Y-junction	900034314	
A	509	Extractor clamp	29250000	
A	510	Glass pane with attachment	29251000	
A	511	Y-Junction for dust extractor	29260000	
A	512	External air duct for SILENT TC/ TC2/ EC2/ powerCAM EC	29370002	
A	513	90° angled connector Silent	29340007	

*) Wear and tear parts, excluded from the guarantee.