

# SHOFU **Universal** Primer

Strong Bonds without Mechanical Retention



# SHOFU Universal Primer

Easy, Efficient and Universal

SHOFU Universal Primer reliably bonds resins to all dental alloys, zirconia and alumina.

In combination with SHOFU Universal Pre-Opaque, SHOFU Universal Primer creates durable bonds, even in a moist environment, and allows you to apply resins to metal frameworks without any mechanical retention.

## Indications

- Application of light-cured opaque resins to frameworks
- Fabrication and repair of resin-veneered restorations
- Intraoral and extraoral repair of porcelain-fused-to-metal restorations
- Bonding of denture resins to metal bases
- Bonding of resins to orthodontic appliances




SHOFU Universal Primer



SHOFU Universal Pre-Opaque

## Benefits

- Quick and easy handling without any additional equipment
- Colourless, no grey areas at the resin-substrate interface
- Reliable bonding, even without the use of retention beads
- For all dental alloys and zirconia-based frameworks
- For intraoral and extraoral use



SHOFU Universal Primer, the successor to the proven SHOFU MZ Primer Plus, is the universal solution for bonding resins to CAD/CAM metal and zirconia frameworks, thanks to its wide range of indications and quick and easy handling without any further accessories or heat treatments.



SHOFU Universal Opaque

**SHOFU Universal Pre-Opaque and SHOFU Universal Primer form a well-matched bonding system, which is extremely resistant to hydrolysis and thermal stress.**

After bonding, simply apply SHOFU Universal Opaque, available in 31 shades, in the tooth shade desired.

# Order Information



## SHOFU Universal Primer

Contents 5 ml, PN Y0060

## Related Products



## SHOFU Universal Pre-Opague

Universal Pre-Opague is highly flowable and bonds strongly to alloys and Universal Opague.

Contents: 2 ml, PN 2111



## SHOFU Universal Opague

Light-cured, ready-to-use opaque paste system, including 16 common tooth shades plus Whitening, NCC and Effect shades.



**SHOFU INC.**, 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan  
**EC REP SHOFU DENTAL GmbH** Am Brühl 17, 40878 Ratingen, Germany, [www.shofu.de](http://www.shofu.de)

SHOFU is a registered trademark of SHOFU INC. All other trademarks and registered trademarks are the property of their respective holders.  
SHOFU INC. reserves the right to change specifications without notice.